



## **Product Specification**

cdhfinechemical.com

	BENEDICT'S REAGENT QUALITATIVE FOR DETECTION OF SUGAR IN URINE		
PRODUCT CODE	808810		
SYNONYMS			
C.I. NO.			
CASR NO.			
ATOMIC ORMOLECULAR FORMULA	-		
ATOMIC ORMOLECULAR WEIGHT	-		
PROPERTIES	_		
PARAMETER	LIMIT		
Description	A clear blue coloured liquid.		
Solubility	Miscible with water.		
Sensitivity test	Passes the test		
	Negative Ative Glucose : Blue 0.5% " : Green 1.0% " : Yellow 1.5% " : Brownish yellow Over 2.0% " : Brick red		
For detection of sugar in urine : Place 5 ml of Bener	dict's qualitative reagent in a test tube. Add 8 drops (0.5 ml) of u	rine. Boil over a flame for 3	
Consult your doctor if you are treated with Strepto Benedict's test for sugar in urine. This solution is not intended for quantitative testin, Note(s) : Assay (if applicable) method menti	g.		
HAZARD WARNING		IMDG Code : 9/III	
Harmful if swallowed. Causes skin irritation. May c	ause an allergic skin reaction. Causes serious eye irritation. May	UN No. : 3082	
cause respiratory irritation. Toxic to aquatic life wit PRECAUTIONARY STATEMENTS		IATA : 9	
	/vanaura / annov Mash _ therewally after handling Do not		
	/ vapours / spray. Wash thoroughly after handling. Do not		
	only outdoors or in a well- ventilated area. Contaminated work ce. Avoid release to the environment. Wear protective		
gloves/protective clothing/eye protection/face pro			
	R or doctor/physician if you feel unwell. IF ON SKIN: Wash with		
	rictim to fresh air and keep at rest in a positioncomfortable fo		
	for several minutes. Remove contact lenses, if present and easy		
to do. Continue rinsing. Call a POISON CENTER or d		y	
unwell. Rinse mouth. If skin irritation occurs: Get n			
	ion. If eye irritation persists: Get medical advice/attention. Take		
	. Wash contaminated clothing before reuse. Collect spillage.		
	e dissolved in dilute acid and transferred to heavy metal waste	]	
Hazard Pictogram(s) :		1	

Replace Date 21-Dec.-2018