

## Technical Information

### Charcoal Blood Agar Base

#### Product Code: DM 1646

**Application:** - Charcoal Blood Agar Base is recommended for the cultivation of *Bordetella pertussis* for vaccine production and also for the maintenance of stock cultures.

#### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Beef extract	10.000
Starch, soluble	10.000
Sodium chloride	5.000
Charcoal	4.000
Yeast extract	3.500
Agar	12.000
Final pH ( at 25°C)	7.5±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

The genus *Bordetella* contains has species : *Bordetella pertussis*, *Bordetella parapertussis*, *Bordetella bronchiseptica* and *Bordetella avium*<sup>(1)</sup>. Genetic studies have shown that these organisms are very closely related to each other. Humans are the only host for *B.pertussis* and *B.parapertussis*, while *B.bronchiseptica* is found in a wide variety of animals and in humans<sup>(2)</sup>. *B.avium* is found in birds. *Bordetella* species are obligately aerobic and metabolically not very active. Except *B.bronchiseptica* all are none-motile .*B.pertussis* is the major cause of whooping cough or pertussis in human beings.

*B.parapertussis* is associated with a milder form of the disease<sup>(3)</sup>. Primary isolation of *B.pertussis* requires an enrichment media in which the addition of charcoal and 15-20% blood neutralize the growth-inhibiting effects of the media.

Charcoal Agar is prepared as per method described by Mishulow, Sharpe and Cohen<sup>(2)</sup>. This medium can be used as a replacement for Bordet-Gengou Agar for isolation of *B.pertussis* from clinical specimens and for the production of *B. pertussis* vaccines. Charcoal Agar supplemented with horse blood can also be used for the cultivation and isolation of *Haemophilus influenzae*<sup>(4)</sup>

The difficulty in the isolation of *Bordetella pertussis* from nasopharyngeal secretions is the inhibition of normal microflora during the long incubation period on nutritious media. Penicillin is added to the medium as an antimicrobial agent for restricting the growth of other contaminants. However Penicillin resistant floras still cause contamination that was observed by Lacey<sup>(4)</sup>. To overcome this problem he supplemented penicillin with diamidino-diphenylamine dihydrochloride, The selectivity of the medium was further increased by adding methicilin by Broome et.al. <sup>(5)</sup>. Sutcliffe and Abbott found that Cephalexin was still better than Methicillin (6) in supressing unwanted nasopharynheal flora associated with clinical specimen.

Regan and Lowe (7) have further showed that Charcoal Blood Agar Base with half strength of cephalexin is an excellent selective enrichment transport medium. Cephalexin is added to inhibit gram-positive organisms that may be present in specimen as cant aminant. Both non-selective and selective media should be inoculated since some stains of *B. pertussis* may be slightly inhibited by cephalexin. Charcoal Blood Agar Base is used for the cultivation of *B.pertussis* for vaccine production.

Medium ingredients like peptic digest of animal tissue, beef extract and yeast extract provide essential nutrients to the organisms. Sodium chloride maintains osmotic balance. Starch soluble and charcoal neutralizes substances toxic to *Bordetella* species such as fatty acids. As charcoal has the tendency to settle at the bottom of the flask swirl the flasks gently to obtain a uniform charcoal suspension while ponring media in petriplates/testtubes<sup>(8)</sup>.

The medium can also be used for the maintenance of stock cultures of *Bordetella pertussis* on slants with weekly subcultures.



Dehydrated Culture Media  
Bases / Media Supplements

## Methodology

Suspend 54.5 grams of powder media in 900 ml distilled water. Shake well & heat to boil to dissolve the medium with frequent stirring. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. For 100 ml of the medium add 10 ml of sterile defibrinated horse blood, 0.3 ml of sterile 100 u/ml Penicillin solution and 0.3 ml of 0.1% solution of 4:4 Diamido-diphenylamine

## Quality Control

### Physical Appearance

Grey to greyish black homogeneous free flowing powder

### Gelling

Firm, comparable with 1.2% Agar gel

### Colour and Clarity of prepared medium

Black coloured, opaque gel with undissolved black particles forms in Petri plates

### Reaction

Reaction of 5.45% w/v aqueous solution at 25°C. pH : 7.5±0.2

pH range 7.30-7.70

### Cultural Response/ characteristics

DM 1646: Cultural characteristics observed w/added sterile defibrinated blood and 100u/ml penicillin solution and 0.1% solution of 4:4 Diamido-diphenylamine hydrochloride, after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
Bordetella bronchiseptica ATCC 4617	50-100	good-luxuriant	>=50%
Bordetella parapertussis ATCC 15311	50-100	good-luxuriant	>=50%
Bordetella pertussis ATCC 8467	50-100	good-luxuriant	>=50%
Staphylococcus aureus ATCC 25923	50-100	inhibited	0%
Klebsiella pneumoniae ATCC 13883	50-100	inhibited	0%

## Storage and Shelf Life

**Dried Media** : Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media** : 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Murray P. R., Baron J. H., Pfaller M. A., Tenover J. C. and Tenover F. C., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Mishulow, Sharpe and Cohen, 1953, Am. J. Public Health, 43:1466
3. Linneman and Pery, 1977, Am. J. Dis. Child., 131:560.
4. Lacey B. W., 1954, J. Hyg., 59:273.
5. Broome C. V., Fraser D. W. and English J. W., 1979, Internat. Symp. on Pertussis DHEW J., Washington D.C., pp 19-29.
6. Sutcliffe E. M. and Abbott J. D., 1979, B.M.J. II: 732-733.
7. Regan and Lowe F., 1977, J. Clin. Microbiol., 6:303.
8. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification -Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Donot use the products if it fails to meet specificatons for identity and performens parameters.





Dehydrated Culture Media  
Bases / Media Supplements



---

**Central Drug House (P) Ltd.** Corp. Office : 7/28, Vardaan House, Ansari Road, Daryaganj, New Delhi-110002 (INDIA)  
Phone : +91-11- 49404040 (100 Lines) Fax : +91-11- 49404050 E-mail : sales@cdhfinechemical.com Website : www.cdhfinechemical.com