



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Basic Fuchsin (6 gm per vial)

Product Code: MS2059

Application: This supplement is recommended for preparation of Endo Agar Base used for confirmation of the presumptive test for the members of the coliform group.

Composition**

Ingredients	Concentration
Basic fuchsin	6g

Directions

Dissolve the contents of 1 vial aseptically in 60 ml of 95% alcohol. Add 4 ml of this (10%) solution to 1000 ml of Endo Agar Base DM2077 / Endo MiVeg Agar Base VM2077 . Heat to boiling to dissolve the medium completely and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile petri plates. Caution: Basic fuchsin is a potential carcinogen and care should be taken to avoid inhalation of the powdered dye and contamination of the skin.

Storage and Shelf Life

Store between 10-30°C and away from light to avoid photo oxidation. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.