

Technical Information

Park and Sanders Selective Supplement II

Product Code: MS2105

Application: An antibiotic supplement recommended for the selective isolation and enumeration of thermotolerant Campylobacter species .

Composition**

Per vial sufficient for 950 ml medium

***Ingredients** Cefoperazone Cycloheximide Concentration 32mg 100mg

Directions

Rehydrate the contents of one vial aseptically with 5 ml of 50% aqueous ethanol. Mix well and aseptically add it to 950 ml of sterile, cooled (45-50ºC) Park and Sanders Broth Base DM2185, DM2185I / Park and Sanders Enrichment MiVeg Broth Base VM2185 alongwith 50 ml of sterile defibrinated lysed horse blood and rehydrated contents of 1 vial of Park and Sanders Selective Supplement I MS2104 / Campylobacter Cefex Agar Base DM2267 containing 50 ml sterile defibrinated sheep blood. Mix well and dispense.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.