



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

M-CP Selective Supplement - II

Product Code: MS2154

Application: Filter sterilized 0.5% solution of phenolphthalein diphosphate, recommended by the Directive of the Council of the European Union 98/83/EC for the selective isolation of *Clostridium perfringens* .

Composition**

Per vial sufficient for 500 ml medium

Ingredients	Concentration
Phenolphthalein diphosphate	0.050g
Distilled water	10ml

Directions

Warm up the refrigerated 0.5% phenolphthalein diphosphate solution to room temperature and add aseptically 10 ml of solution to 485 ml sterile, molten, cooled (45-50°C) M-CP Agar Base DM2354 / M-CP MiVeg Agar Base VM2354 along with rehydrated contents of one vial of M-CP Selective Supplement I MS2153 . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.