

# **Technical Information**

### **Cefoxitin Supplement**

Product Code: MS2259

**Application:** An antimicrobial supplement recommended for the selective isolation of Methicillin Resistant *Staphylococcus aureus* from clinical specimens.

### Composition\*\*

Per vial sufficient for 500 ml medium

\*Ingredients Concentration
Cefoxitin 3mg

#### **Directions**

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add to 500 ml of sterile molten MeReSa Agar Base DM2594 / MiCrome MeReSa Agar, Base DM2674 / MiCrome MeReSa MiVeg Agar Base VM1674 / MiCrome MRSA Agar Base, Modified DM2953 . This supplement can either be used individually or in combination with MS2229 MeReSa Selective Supplement for more selectivity. Mix well and pour into sterile petri plates

## Storage and Shelf Life

Store at 2-8°C. Use before the expiry on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.