



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Staph-Strepto Supplement

Product Code: MS2030

Application: An antibiotic supplement recommended for the selective isolation of Staphylococci and Streptococci species.

Composition**

Per vial sufficient for 500 ml medium

| *Ingredients | Concentration |
|-------------------|---------------|
| Nalidixic acid | 1mg |
| Colistin sulphate | 1mg |

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base DM1144 or Columbia Blood Agar Base, MiVeg VM1144 along with 25 ml of sterile defibrinated horse blood. Mix gently and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.