



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Mycoplasma Enrichment Supplement

Product Code: MS2075

Application: A selective supplement for the isolation of Mycoplasma.

Composition**

Per vial sufficient for 70 ml medium

*Ingredients	Concentration
Horse serum	20ml
Yeast extract (25%w/v solution)	10ml
Thallos acetate	25mg
Penicillin G	20000IU

Directions

Aseptically add the contents of 1 vial to 70 ml of sterile, molten, cooled (45-50°C) Mycoplasma Agar Base DM1266 / Mycoplasma MiVeg Agar Base VM1266 / Mycoplasma Broth Base w/o CV DM1267 / Mycoplasma MiVeg Broth Base w/o CV VM1267 / Mycoplasma Broth Base w/ CV DM1268 / Mycoplasma MiVeg Broth Base w/ CV VM1268 / PPLO Modified Broth Base w/o CV DM2586 . For Mycoplasma Broth Base w/ CV DM1268 / Mycoplasma MiVeg Broth Base w/ CV VM1268 add 0.285 ml of 1% Potassium Tellurite MS2052 along with Mycoplasma Enrichment Supplement MS2075 . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.