

# **Technical Information**

### \*\*L.mono Enrichment Supplement I

## Product Code: MS2214

Application: For selective differentiation of Listeria monocytogenes from other Listeria species, as per ISO Committee.

#### Composition\*\*

Per vial sufficient for 500 ml medium

Ingredients	Concentration
L – phosphatidylinositol	1 g
Distilled water	25 ml

#### Directions

Thaw the contents of 1 vial of L. mono Enrichment Supplement I at room temperature. Aseptically add the sterile contents to 460 ml of sterile molten, cooled (45-50°C) L. mono Differential Agar Base DM2540 / L. mono Differential MiVeg Agar Base VM2540 along with sterile rehydrated contents of 1 vial each of L. mono Selective Supplement I MS2212 and L. mono Selective Supplement II MS2213 . Mix well and pour into sterile petri plates.

### Storage and Shelf Life

On receipt product should be stored at -20°C. Use before the expiry date on the label.

### **Disclaimer**

- User must ensure suitability of the product(s) in their application prior to use. •
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time. •
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for • infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters. ٠