

Technical Information

Microme Selective Salmonella Agar Supplement

Product Code: MS2274

Application: For the selective isolation and differentiation of *Salmonella* species from coliforms by chromogenic method.

Composition**

Per vial sufficient for 1000 ml medium

*Ingredients	Concentration
Novobiocin	10mg
Cefsulodin	24mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 1000 ml of sterile, cooled MiCrome Selective Salmonella Agar Base DM2842. Mix well and pour into sterile Petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.