

## Technical Information

### Enterococcus faecium Selective Supplement

**Product Code: MS2226**

**Application:** Recommended to differentiate *Enterococcus faecium* from *Enterococcus faecalis*.

#### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Cephalexin	25mg
Aztreonam	37.500mg

#### Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Arabinose Agar Base DM2576 / MiCrome Enterococcus faecium Agar Base DM2580 / MiCrome Enterococcus Faecium MiVeg Agar Base VM2580 . Mix well and pour into sterile petri plates.

#### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.