



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

T.S.C. Supplement (Perfringens T.S.C. Supplement)

Product Code: MS2014

Application: An antibiotic supplement, recommended for the selective isolation of *Clostridium perfringens*.

Composition**

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
D-Cycloserine	200mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add to 475 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base (T.S.C. / S.F.P., DM1837), Perfringens HiVeg™ Agar Base (T.S.C. / S.F.P., VM1837), alongwith 25 ml of Egg Yolk Emulsion MS2045 or 500 ml of Perfringens Agar Base DM8371 / Tryptose Cycloserine Dextrose Agar Base DM2233 / Tryptose Cycloserine Dextrose MiVeg Agar Base VM2233 / Tryptose Cycloserine Azide Agar Base DM2279 . Tryptone Yeast Sodium sulphite Agar Base DM2046I.

Mix well and pour into sterile plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.