



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Basic Fuchsin (0.17 gm per vial)

Product Code: MS2059A

Application: Recommended for preparation of Endo Agar which is used for confirmation of the presumptive test for the members of the coliform group.

Composition**

Per vial sufficient for 1000 ml medium

Ingredients	Concentration
Basic fuchsin	0.170g

Directions

Dissolve contents of 1 vial in 10 ml of 95% ethanol. After trituration for dissolving the dye, filter the solution if necessary. Add this solution to 990 ml of Modified Endo Agar Base DM2481. Heat to boiling to dissolve the medium completely and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates. Caution : Basic fuchsin is a potential carcinogen and care should be taken to avoid inhalation of powdered dye and contamination of the skin.

Storage and Shelf Life

Store between 10-30°C and away from light to avoid photo oxidation. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.