

Technical Information

Tributyrin (10 ml per vial)

Product Code: MS2081

Application: Recommended for detection of lipolytic microorganisms.

Composition**

Per vial sufficient for 990 ml medium
Ingredients Concentration
Tributyrin 10ml

Directions

Warm up the refrigerated solution to room temperature. Add contents of one vial aseptically to 990 ml of Tributyrin Agar Base w/o Tributyrin DM1157 / Tributyrin MiVeg Agar Base w/o Tributyrin MV1157 and then sterilize by autoclaving at 15 lbs pressure for 15 minutes. Mix well and pour into sterile petri plates

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.