



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Perfringens Supplement-II

**Product Code:MS2012**

**Application:** An antibiotic supplement, recommended for the selective isolation of *Clostridium perfringens* .

### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Oleandomycin phosphate	0.250mg
Polymyxin B sulphate	5000IU

### Directions

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base O.P.S.P. DM1579 / Perfringens MiVeg Agar Base O.P.S.P. VM1579 alongwith sterile rehydrated contents of 1 vial of Perfringens Supplement I MS2011 . Mix well and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.