



Product Specification

cdhfinechemical.com

POTASSIUM PERMANGANATE SOLUTION 0.02 MOL/LIT CITRISOL 1.5805g KMnO₄ FOR 500 ML SOLUTION 0.1N SOLUTION

PRODUCT CODE	878780
SYNONYMS	
C.I. NO.	-
CASR NO.	-
ATOMIC OR MOLECULAR FORMULA	
ATOMIC OR MOLECULAR WEIGHT	
PROPERTIES	
PARAMETER	LIMIT
Description	A clear purple coloured solution.
Solubility	Miscible with water.
Normality	0.1N ± 0.002
(After diluting total solution to 500 ml with distilled water, at 27°C.)	

Direction for use and storage: Ampoules

- 1. Hold this glass ampoule firmly in hand to prevent it from falling while cutting the tip. With a sharp knife, slice offthe tip. This will break the tip smoothly the rim.
- 2. Transfer the contents of the ampoule into a 500 ml standard volumetric flask of Class A grade. Rinse the ampouleand the tip carefully with distilled water. Transfer the rinsings to the flask. A small funnel provided with pack may beused for this purpose. Make up volume up to the mark with distilled water at 27°C. Stopper the flask and mixthoroughly.
- 3. Note down the normality and date of preparation of the solution on the sticker label provided in the pack. Affix thesame on the bottle / flask.
- 4. The unopened ampoules may be stored indefinitely without the contents deteriorating in strength or quality.(Solutions diluted for use should be kept under conditions of storage appropriate to volumetric and reagentsolutions.)
- 5. The accuracy of the solution prepared as above from SD's Concentrated Volumetric Solution is 0.1 ± 0.002 N.

Note(s): Assay (if applicable) method mentioned.

WARNING
Hazard statements: Causes mild skin irritation. Causes eye irritation. Toxic to aquatic life with long lasting effects.

Precautionary statements:

Prevention: Wash hands thoroughly after handling.

IMDG Code: 9/III
UN No.: 3082
IATA: 9

Response: If skin irritation occurs, seek medical advice/ attention. If eye irritation persists, get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, If present and

advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, If present and easy to do. Continue rinsing.

Disposal: Dispose of contents and container in accordance with relevant legislation.

Hazard Pictogram(s):



Replace Date 26-Nov.-2018