

Dehydrated Culture Media Bases / Media Supplements

# **Technical Information**

### **Meat Extract Powder**

### Product Code: BA2003

### **Principle & Interpretation**

Meat Extract Powder is specially selected to yield maximum growth of fastidious microorganisms. It is a yellowish brown coloured powder having an agreeable meat like odour and taste. It readily dissolves in distilled water, forming yellowish brown coloured solution, which remains clear even after autoclaving. It is recommended for general purpose and diagnostic media preparations.

### Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Safety guidelines may be referred in individual safety data sheets.

#### Limitations

1. It is biological origin product since variation in colour of powder and clarity may observed.

2. Each lot of the product has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user 's requirement.

3. Individual organisms differ in their growth requirement and may show variable growth patterns on themedium prepared by the product.

#### Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

## **Quality Control**

#### Appearance

Light yellow to brownish yellow homogenous free flowing powder characteristic odour but not putrescent

#### Solubility

Freely soluble in distilled/purified water, insoluble in alcohol.

#### Clarity

1% w/v aqueous solution remains clear without haziness after autoclaving at 15 lbs pressure (121ºC) for15 minutes.

#### рΗ

pH of 2% w/v aqueous solution at 25ºC 5.9-6.9

#### Microbial Load

Bacterial Count : <= 2000 CFU/gram by plate method, when incubated at 30-35ºC for not less than 3 days Yeast & mould Count : <= 100 CFU/gram by plate method, when incubated at 20-25ºC for not less than 5 days.

#### Test for pathogens

- 1. Escherichia coli- Absent/gram of sample
- 2. Salmonella species- Absent/10 gram of sample
- 3.Pseudomonas aeruginosa- Absent/gram of sample
- 4. Staphylococcus aureus- Absent/gram of sample
- 5. Candidaalbicans- Absent/gram of sample
- 6. Clostridia- Absent/gram of sample
- Indole Test : Tryptophan content: Passes

Cultural response : Cultural response observed after an incubation at 35-37°C for 18-24 hours by preparing Sterility Medium A



Dehydrated Culture Media Bases / Media Supplements

(DM1017) using Extract Powder as an ingredient. Cultural response Organism

OrganismGrowthBacillus subtilis subsp. Spizizenii ATCC 6633(WDCM 00003)LuxuriantPseudomonas aeruginosa ATCC 27853(WDCM 00025)LuxuriantStaphylococcus aureus subsp.aureus ATCC25923(WDCM 00034)LuxuriantEnterococcus faecalis ATCC 29212 (WDCM 00087)Luxuriant

#### Chemical Analysis :

Total Nitrogen : ≥12.50 %

Amino Nitrogen : ≥3.50 %

Sodium chloride : ≤5.50 %

Loss on drying : ≤6.00 %

Residue on ignition : ≤12.00 %

### Storage and Shelf Life

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

#### Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

### **Disclaimer :**

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
  Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.