

Technical Information

Brain Heart Infusion Broth

Product Code: DM 1210

Application: Brain Heart Infusion Broth is employed for the propagation of fastidious pathogenic cocci and other organisms associated with blood culture work and allied pathological investigations.

Composition**

| Ingredients | Gms / Litre |
|---------------------------|-------------|
| HM infusion powder # | 12.500 |
| Beef heart, infusion from | 5.000 |
| Proteose peptone | 10.000 |
| Dextrose | 2.000 |
| Sodium chloride | 5.000 |
| Disodium phosphate | 2.500 |
| Final pH (at 25°C) | 7.4±0.2 |

**Formula adjusted, standardized to suit performance parameters

Equivalent to Calf brain infusion from

Principle & Interpretation

Brain Heart Infusion Medium is useful for cultivating a wide variety of microorganisms and preparing the inocula for antimicrobial susceptibility testing. This medium is a modification of the original formulation of Rosenow, where he added pieces of brain tissues to dextrose broth ⁽¹⁾. Brain Heart Infusion Broth is also preferred for growing anaerobic bacteria, yeasts and moulds ⁽²⁻⁴⁾. This medium is nutritious and well buffered to support the growth of wide variety of organisms ^(2, 5, 6). Addition of 10% of defibrinated sheep blood, make this media useful for isolation and cultivation of *Histoplasma capsulatum* ⁽⁷⁾. Also addition of gentamicin and/or chloramphenicol is recommended for isolation of other fungi ⁽⁸⁾. Proteose peptone and infusions (calf brain and beef heart) serve as sources of carbon, nitrogen, essential growth factors, amino acids and vitamins. Dextrose serves as a source of energy. Disodium phosphate helps in maintaining the buffering action of the medium whereas sodium chloride maintains the osmotic equilibrium of the medium.

Methodology

Suspend 37 grams of powder media in 1000 ml distilled water. Shake well & dispense into bottles or tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. For best results, the medium should be used on the day it is prepared, otherwise, it should be boiled or steamed for a few minutes and then cooled before use.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light to medium amber coloured, clear solution without any precipitate

Reaction

Reaction of 3.7% w/v aqueous solution at 25°C. pH : 7.4±0.2

pH range 7.20-7.60

Cultural Response/Characteristics

DM 1210: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

| Organism | Inoculum (CFU) | Growth |
|--|----------------|----------------|
| <i>Enterococcus faecalis</i> ATCC 29212 | 50-100 | good-luxuriant |
| <i>Neisseria meningitides</i> ATCC13090 | 50-100 | good-luxuriant |
| <i>Streptococcus pneumoniae</i> ATCC 6303 | 50-100 | good-luxuriant |
| <i>Streptococcus pyogenes</i> ATCC 19615 | 50-100 | good-luxuriant |
| <i>Candida albicans</i> ATCC10231 | 50-100 | good-luxuriant |
| <i>Staphylococcus aureus</i> ATCC 25923 | 50-100 | good-luxuriant |

Storage and Shelf Life

Dried media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Rosenow, 1919, J. Dental Research, 1:205.
2. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.
3. Atlas R. M., 1993, Handbook of Microbiological Media, 147-153, CRC Press, Boca Raton, FL.
4. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
5. Roseburg T. et al, 1944, J. Inf. Dis., 74:13 1
6. Conant N. F., 1950, Diagnostic Procedures and Reagents, 3rd Ed., APHA Inc., New York
7. Howard B., Keiser J. F., Weissfeld A. et al, 1994, Clinical and Pathogenic Microbiology, 2nd Ed., Mosby Co.
8. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

Replaced Date 08-Jan-2026