

Dehydrated Culture Media Bases / Media Supplements

Technical Information

Special Peptone, Certified

Product Code: BC2015

Principle & Interpretation

Special Peptone, Certified is designed to maximize the growth of fastidious microorganisms. It is carefully processed to ensure consistent growth characteristics when used in differential media.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. It is biological origin product since variation in colour of powder and clarity may observed.

2. Each lot of the product has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's requirement.

3. Individual organisms differ in their growth requirement and may show variable growth patterns on themedium prepared by the product.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

Quality Control

Appearance

Yellow to brownish yellow homogenous free flowing powder characteristic but not putrescent

Solubility

Freely soluble in distilled/purified water, insoluble in alcohol.

Clarity

1% w/v aqueous solution is clear without any haziness after autoclaving at 15 lbs pressure (121°C) for 15minutes.

рΗ

pH of 2% w/v aqueous solution at 25⁰C Microbial Load 6.1 - 7.1

Bacterial Count : <= 2000 CFU/gram by plate method, when incubated at 30-35⁰C for not less than 3 days Yeast & mould Count : <=

100 CFU/gram by plate method, when incubated at 20-25ºC for not less than 5 days.

Test for pathogens

- 1. Escherichia coli- Absent/gram of sample
- 2. Salmonella species- Absent/10 gram of sample
- 3. Pseudomonas aeruginosa- Absent/gram of sample
- 4. Staphylococcus aureus- Absent/gram of sample
- 5. Candidaalbicans- Absent/gram of sample
- 6. Clostridia- Absent/gram of sample

Indole Test : Tryptophan content: Passes



Dehydrated Culture Media Bases / Media Supplements

Cultural response : Cultural response observed after an incubation at 35-37°C for 40-48 hours by preparingColumbia Agar Base (DM1144) using Special Peptone, Certified as an ingredient.

Cultural response			
	Organism	Growth	Haemolysis
	Neisseria meningitidis ATCC 13090	Luxuriant	none
	Staphylococcus aureus subsp. aureus ATCC25923 (WDCM 00034)	Luxuriant	beta / gamma
	Staphylococcus epidermidis ATCC 12228(WDCM 00036)	Luxuriant	gamma
	Streptococcus pneumoniae ATCC 6303	Luxuriant	alpha
	Streptococcus pyogenes ATCC 19615	Luxuriant	beta

Chemical Analysis

Total Nitrogen : ≥11.50 %

Amino Nitrogen : ≥3.00 %

Sodium chloride : ≤7.00 %

Loss on drying : ≤5.00 %

Residue on ignition : ≤15.00 %

Storage and Shelf Life

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use .

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing
 of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.