

# **Technical Information**

## **Meat Peptone Type P**

Product Code: BA3049

# Principle & Interpretation

HM Peptone Type P is processed under controlled conditions to retain all nutritive values. It is obtained by proteolytic activity of Pepsin. It yields heavy growth of fastidious microorganisms when used in nutrient and differential media. Also recommended for industrial fermentation studies and production of antibiotics, secondary metabolites, enzymes etc.

## Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Safety guidelines may be referred in individual safety data sheets.

#### Limitations

- 1. It is biological origin product since variation in colour of powder and clarity may observed.
- 2. Each lot of the product has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's requirement.
- 3. Individual organisms differ in their growth requirement and may show variable growth patterns on themedium prepared by the product.

#### Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

# **Quality Control**

#### Appearance :

Light yellow to brownish yellow homogenous free flowing characteristic odour but not putrescent.

#### Solubility:

Freely soluble in distilled/purified water, insoluble in alcohol and ether.

#### Clarity:

1% w/v aqueous solution is clear without any haziness after autoclaving at 15 lbs pressure (121ºC) for 15minutes.

#### pH:

pH of 2% w/v aqueous solution at 25°C=6.0 - 7.0

#### Microbial Load:

Bacterial Count : <= 2000 CFU/gram by plate method, when incubated at 30-35°C for not less than 3 days. Yeast & mould Count : <= 100 CFU/gram by plate method, when incubated at 20-25°C for not less than 5 days.

### Test for pathogens:

- 1. Escherichia coli- Absent/gram of sample
- Salmonella species- Absent/10 gram of sample
- 3. Pseudomonas aeruginosa- Absent/gram of sample
- 4. Staphylococcus aureus- Absent/gram of sample
- 5. Candidaalbicans- Absent/gram of sample
- 6. Clostridia- Absent/gram of sample

Indole Test: Tryptophan content: Passes

Cultural response: Cultural response observed after an incubation at 35-37°C for 18-24 hours byreparing Nutrient Agar



(DM1001), using HM Peptone Type P as an ingredient.

### **Cultural response**

Organism	Growth
Escherichia coli ATCC 25922 (WDCM00013)	Luxuriant
Pseudomonas aeruginosa ATCC 27853(WDCM 00025)	Luxuriant
Staphylococcus aureus subsp.aureus ATCC25923 (WDCM 00034)	Luxuriant
Salmonella enterica subsp. enterica TyphiATCC 6539	Luxuriant
Streptococcus pyogenes ATCC 19615	Luxuriant
Salmonella enterica subsp.enterica EnteritidisATCC 13076 (WDCM 00030)	Luxuriant
Salmonella enterica subsp.enterica Typhimurium ATCC 14028 (WDCM 00031)	Luxuriant
Yersinia enterocolitica ATCC 9610(WDCM 00038)	Luxuriant
Yersinia enterocolitica ATCC 23715(WDCM 00160)	Luxuriant

#### Chemical Analysis

Total nitrogen: ≥13.00 %

Amino Nitrogen: ≥3.50 %

Sodium Chloride: ≤6.00 %

Loss on drying: ≤5.00 %

Residue on ignition: ≤12.00 %

### Storage and Shelf Life

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

# Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

### Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.