

# **Technical Information**

## Peptone Type I, Bacteriological

Product Code: BA2667

## Principle & Interpretation

Peptone Type I, Bacteriological is prepared under controlled conditions, used as culture media ingredient for cultivation of various bacteria. Employed for different fermentation processes, production of vaccines, antibiotics, etc.

### Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Safety guidelines may be referred in individual safety data sheets.

#### Limitations

- 1. It is biological origin product since variation in colour of powder and clarity may observed.
- 2. Each lot of the product has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user 's requirement.
- 3. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium prepared by the product.

#### Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature

## **Quality Control**

#### Appearance

Light yellow to brownish yellow homogeneous free flowing powder characteristic odour but not putrescent

#### Solubility

Freely soluble in distilled/purified water, insoluble in alcohol.

#### Clarity

1% w/v aqueous solution is clear to slight hazy may develop precipitate after autoclaving at 15 lbs pressure(121ºC) for 15 minutes.

#### На

pH of 2% w/v agueous solution at 25°C. 6.50-7.50

#### Microbial Load :

Bacterial Count: <= 10000 CFU/gram by plate method, when incubated at 30-35°C for not less than 3 days Yeast & mould Count: <= 100 CFU/gram by plate method, when incubated at 20-25°C for not less than 5 days.

#### Test for pathogens:

- 1. Escherichia Coli Absent/gram of sample
- 2. Salmonella species- Absent/10 gram of sample
- 3. Pseudomonas aeruginosa Absent/gram of sample
- 4. Staphylococcus aureus Absent/gram of sample
- 5. Candidaalbicans Absent/gram of sample



6. Clostridia- Absent/gram of sample

Indole Test: Tryptophan content: Passes

Cultural response: Cultural response observed after an incubation at 35-37°C for 18-24 hours by preparingNutrient Agar (DM1001) using Peptone Type I, Bacteriological as an ingredient.

#### Cultural response

Organism	Growth
Escherichia coli ATCC 25922 (WDCM00013)	Luxuriant
Pseudomonas aeruginosa ATCC 27853(WDCM 00025)	Luxuriant
Staphylococcus aureus subsp. aureus ATCC25923 (WDCM 00034)	Luxuriant
Salmonella enterica subsp. enterica TyphiATCC 6539	Luxuriant
Streptococcus pyogenes ATCC 19615	Luxuriant
Salmonella enterica subsp. Enterica Enteritidis ATCC 13076 (WDCM 00030)	Luxuriant
Salmonella enterica subsp. Enterica Typhimurium ATCC 14028(WDCM 00031)	Luxuriant
Yersinia enterocolitica subsp. Enterocolitica ATCC 9610 (WDCM 00038)	Luxuriant
Yersinia enterocolitica subsp. Enterocolitica ATCC 23715 (WDCM 00160)	Luxuriant

#### Chemical Analysis

Total Nitrogen : ≥11.30 %

Amino Nitrogen : ≥3.00 %

Sodium chloride : ≤6.00 %

Loss on drying : ≤5.00 %

Residue on ignition : ≤15.00 %

### Storage and Shelf Life

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use

### Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

### Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.