

## Technical Information

### Violet Red Bile Glucose Agar

**Product Code: DM 1581 H**

**Application:** Violet Red Bile Glucose Agar is recommended for detection and enumeration of *Enterobacteriaceae* from pharmaceutical products in accordance with the microbial limit testing by harmonized methodology of USP/EP/BP/JP/IP (Medium 5).

### Composition\*\*

Ingredients	Gms / Litre
Yeast extract	3.000
Pancreatic digest of gelatin	7.000
Bile salts	1.500
Sodium chloride	5.000
Glucose monohydrate	10.000
Agar	15.000
Neutral red	0.030
Crystal violet	0.002
pH after heating ( at 25°C)	7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Principle & Interpretation

Violet Red Bile Glucose Agar is a selective medium recommended for detection and enumeration of members of *Enterobacteriaceae* especially the bile tolerant gram negative bacteria in accordance with the microbial limit testing and as per the methodology of USP/EP/BP/JP/IP (1,2,3,4,5) from non-sterile products and pharmaceutical preparations.

Pancreatic digest of gelatin and yeast extract provide nitrogenous compounds and other nutrients essential for bacterial metabolism. Bile salts and crystal violet make this media selective. Crystal violet inhibits gram-positive organisms especially Staphylococci. Neutral red indicator helps to detect glucose fermentation. Glucose fermenting strains produce red colonies with pink-red halos in the presence of neutral red. Sodium chloride maintains the osmotic equilibrium in the medium. The red colour is due to absorption of neutral red and a subsequent colour change of the dye when the pH of medium falls below 6.8.

### Methodology

Suspend 40.62 grams of powder media (the equivalent weight of dehydrated medium per litre) in 1000 ml purified /distilled water. Shake well & heat to dissolve the medium completely. DO NOT HEAT IN AN AUTOCLAVE. Mix well and pour into sterile Petri plates.

### Quality Control

#### Physical Appearance

Light yellow to pinkish beige homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.5% Agar gel.

#### Colour and Clarity of prepared medium

Reddish purple coloured clear to slightly opalescent gel forms in Petri plates.

#### pH Range:-

7.20-7.60

#### Growth Promotion Test

Growth Promotion is carried out in accordance with the harmonized method of USP/EP/BP/JP. Cultural response was observed after an incubation at 30-35°C for 18-24 hours. Recovery rate is considered as 100% for bacteria growth on Soyabean Casein Digest Agar.

#### Growth promoting properties

Growth of microorganism comparable to that previously obtained with previously tested and approved lot of medium occurs at the specified temperature for not more than the shortest period of time specified inoculating  $\leq 100$  cfu (at 30-35°C for  $\leq 18$  hours).

### Cultural Response/Characteristics

DM 1581: Colonies are comparable in appearance and indication reaction to those previously obtained with previously tested and approved lot of medium occurs for the specified temperature for a period of time within the range specified inoculating  $\leq 100$  cfu (at 30-35°C for 18-24 hours).

Organism	Inoculum (CFU)	Growth	Observed Lot value (CFU)	Recovery	Colour of colony	Incubation temperature
<i>Escherichia coli</i> ATCC 8739	50-100	luxuriant	25 -100	$\geq 50$ %	pink-red with bile precipitate	18 -24 hrs
<i>Pseudomonas aeruginosa</i> ATCC 9027	50-100	Luxuriant	25 -100	$\geq 50$ %	pink to red	18 -24 hrs
<b>Additional Microbiological Testing</b>						
<i>Escherichia coli</i> NCTC 9002	50-100	good-luxuriant	25 -100	$\geq 50$ %	pink-red with bile precipitate	18 -24 hrs
<i>Escherichia coli</i> ATCC 25922	50-100	good-luxuriant	25 -100	$\geq 50$ %	pink-red with bile precipitate	18 -24 hrs
<i>Salmonella Enteritidis</i> ATCC 13076	50-100	good-luxuriant	25 -100	$\geq 50$ %	light pink	18 -24 hrs
<i>Enterobacter aerogenes</i> ATCC 13048		good-luxuriant	25 -100	$\geq 50$ %	pink-red	18 -24 hrs
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	Inhibited	25 -100	0%		$\geq 24$ hrs
<i>Staphylococcus aureus</i> ATCC 6538	$\geq 10^3$	Inhibited	25 -100	0%		$\geq 24$ hrs

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. The United States Pharmacopoeia, 2011, The United States Pharmacopoeial Convention. Rockville, MD.
2. British Pharmacopoeia, 2011, The Stationery office British Pharmacopoeia 7
3. European Pharmacopoeia, 2011, European Dept. for the quality of Medicines.
4. Japanese Pharmacopoeia, 2008.
5. Indian Pharmacopoeia, 2010 Ministry of Health and Family Welfare, Govt. of India

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

Replace date 03-July-2025