

Technical Information

Salmonella Selective Enrichment Supplement

Product Code: MS2275

Application: Selective supplement recommended for Selective isolation and differentiation of Salmonella species from coliforms by chromogenic method.

Composition**

Per vial sufficient for 1000 ml medium

| *Ingredients | Concentration |
|--------------|---------------|
| Novobiocin | 12mg |

Directions

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents of one vial to 1000ml of sterile, cooled (45-50°C) of Salmonella Selective Enrichment Broth Base DM2843 . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.