

Technical Information

PP Pseudomonas Selective Supplement II

Product Code: MS2265

Application: An antimicrobial supplement recommended for the selective isolation of *Pseudomonas* species.

Composition**

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Pimaricin (natamycin)	0.010g

Directions

Rehydrate the contents of 1 vial with 10 ml of sterile distilled water and aseptically to 1000 ml of sterile molten Pimaricin and Pseudomonas Selective Supplement (PP Pseudomonas Agar Base) DM2788 along with 1 vial of PP Pseudomonas Selective Supplement MS2264 . Mix well and pour into sterile Petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.