

Technical Information

Potato Dextrose Broth, Granulated

Product Code: GM 1403

Application: - Potato Dextrose Broth, granulated is recommended for the isolation and enumeration of yeasts and moulds.

Composition**

Ingredients	Gms / Litre	
Potatoes, infusion from	200.000	
Dextrose	20.000	
Final pH (at 25°C)	5.1±0.2	

^{**}Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Potato Dextrose Broth is recommended by APHA (1) and F.D.A. (2) for plate counts of yeasts and moulds in the examination of foods and dairy products (3). Potato Dextrose Broth is also used for stimulating sporulation, for maintaining stock cultures of certain dermatophytes and for differentiation of typical varieties of dermatophytes on the basis of pigment production (4). Potato infusion and dextrose promote luxuriant fungal growth. Adjusting the pH of the medium by tartaric acid to 3.5, inhibits the bacterial growth. Heating the medium after acidification should be avoided.

Methodology

Suspend 24.0 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well before dispensing. In specific work, when pH 3.5 is required, acidify the medium with sterile 10% tartaric acid. The amount of acid required for 100 ml. of sterile, cooled medium is approximately 1 ml. Do not heat the medium after addition of the acid.

Quality Control

Appearance

Off-white to yellow coloured granular medium.

Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent solution in tubes.

Reaction

Reaction of 2.4% w/v aqueous solution at 25°C. pH: 5.1±0.2

pH Range

4.90-5.30

Cultural Response

Cultural characteristics observed after an incubation at 25-30°C for 4-5 days.

Organism	Inoculum (CFU)	Growth	Ascospore formation
Cultural Response #Aspergillus brasiliensis ATCC 16404 (00053*)	50-100	luxuriant	negative
Candida albicans ATCC 10231 (00054*)	50-100	luxuriant	negative
Saccharomyces cerevisiae ATCC 9763 (00058*)	50-100	luxuriant	positive





Key:-# Formerly known asAspergillus niger.

Key:- * Corresponding WDCM number.

Storage and Shelf Life

Dried Media: On receipt store between 10-30°C in tightly closed container and the prepared medium at 2 - 8°C.

Prepared Media: Use before expiry date on the label.

Further Reading

- 1. Salfinger Y. and Tortorello, ML (Eds.), 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., APHA, Washington, D.C.
- 2. FDA Bacteriological Analytical Manual, 2005, 18th Ed., AOAC, Washington, DC.
- 3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
- 4. MacFaddin J. F., 1985, Media for the Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol.1, Williamsand Wilkins, Baltimore

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for
 infringement of any patents. Do not use the products if it fails to meet specificatons for identity and performens parameters.

