



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Mycoplasma Supplement - P

Product Code: MS2182

Application: A selective supplement for the isolation of *Mycoplasma pneumoniae* when used with Mycoplasma Agar Base / Mycoplasma MiVeg Agar Base.

Composition**

*Ingredients	Concentration
Horse serum	6ml
Yeast extract	3 gm
Thallos acetate	0.008g
Penicillin	12000IU
Glucose	0.300g
Phenol red	0.0012g
Methylene blue	0.0003g
Mycoplasma Broth Base w/o CV (DM1267)	0.145g

Directions

Aseptically add 20 ml sterile distilled water to the vial and mix gently. Separately prepare Mycoplasma Agar Base DM1266 /Mycoplasma MiVeg Agar Base VM1266 without supplement addition. Heat to boiling to dissolve the medium completely.

Add 1 ml of the medium to each 10 small bottles. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. After sterilization allow the medium to set or solidify. Then aseptically add 2 ml of the reconstituted Mycoplasma Supplement-P MS2182 to each bottle containing agar. This will form a biphasic medium.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.