



Ready Prepared Media

Technical Information

Sheep Blood Agar, Modified Plate

Product Code: PM 2956

Application: Recommended for cultivation and studying haemolytic reactions of *Bacillus cereus*. The composition and performance of this media are as per the specification laid down in ISO 21871 :2006

Composition**

Ingredients	Gms / Litre
Tryptone ##	15.000
Soya peptone#	5.000
Sodium chloride	5.000
Agar	12.500
Blood	50ml
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

#- Equivalent to Enzymatic digest of soya

##- Equivalent to Enzymatic digest of casein

Principle & Interpretation

Haemolysins are exotoxins produced by bacteria that lyse red blood cells. The haemolytic reaction can be visualized on blood agar plates. On blood agar plates colonies of haemolytic bacteria may be surrounded by clear, colourless zone where the red blood cells have been lysed and the haemoglobin destroyed to a colourless compound. This is beta haemolysis. Other types of bacteria can reduce haemoglobin to methaemoglobin which produces a greenish zone around the colonies and is called alpha haemolysis (5). Gamma haemolysis is no haemolysis where no change in the medium is observed (3). *Bacillus cereus* is Gram -positive aerobic or facultatively anaerobic, motile, spore forming, rod shaped bacterium that is widely distributed environmentally. *B.cereus* is associated mainly with food poisoning it is increasingly reported to be cause of serious and fatal non- gastrointestinal-tract infections.

Sheep Blood Agar Base, Modified with added sheep blood was developed to allow maximum recovery of *B.cereus* without interfering with their haemolytic reactions. This medium is formulated in accordance with ISO (4). It was formulated to be compatible with sheep blood and give improved haemolytic reactions of organisms. Tryptone and Soya peptone provide nitrogen, carbon, amino acids and vitamins. Sodium chloride maintains the osmotic balance.

Type of specimen

Food and animal feeding stuff samples.

Specimen Collection and Handling

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (4,6).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :



Ready Prepared Media

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium .
2. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Methodology

Either streak, inoculate or surface spread the test inoculum (50-100 CFU) aseptically on the plate.

Quality Control

Appearance

Sterile Sheep Blood Agar, Modified in 90 mm disposable plates.

Colour

Cherry coloured medium

Quantity

25 ml of medium in 90mm disposable plate

pH

7.10-7.50

Sterility Check

Passes release criteria

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours

Organism	Growth	Inoculum(CFU)	Recovery	Haemolysis
Bacillus cereus ATCC 10876	luxuriant	50-100	>=70%	beta

Storage and Shelf Life

- On receipt store between 2-8°C Use before expiry date on the label.
- Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (1,2).



Ready Prepared Media

Further Reading

1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
2. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
3. Koneman E. W., Allen S. D., Janda W. M., Schreckenberger P. C., Winn W. C. Jr., 1992, Colour Atlas and Textbook of Diagnostic Microbiology, 4th Ed., J. B. Lippincott Company.
4. Microbiology of Food & Animal feeding stuffs. Horizontal method for the determination of low numbers of presumptive *Bacillus cereus*- Most probable number technique and detection methods International Organization for Standardization (ISO), Draft ISO 21871:2006 .
5. Pelczar M. J. Jr., Reid R. D., Chan E. C. S., 1977, Microbiology, 4th Ed., Tata McGraw-Hill Publishing Company Ltd, New Delhi.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.