

### Technical Information

#### Phyto Xano Camp Agar Base

#### Product Code: PHM1019

**Application:** Selective medium for the isolation of *Xanthomonas campestris pv campestris* and *Xanthomonas campestris pv. armoraciae* from crucifer seeds.

#### Composition\*\*

Ingredients	Grams/Litre
Soya Peptone	2.00
Tryptone	2.00
Soluble Starch	25.00
Potassium dihydrogen phosphate	2.80
Diammonium hydrogen phosphate	0.80
Magnesium Sulphate anhydrous	0.1952
L-Glutamine	6.00
L- Histidine	1.00
Glucose monohydrate	1.00
Agar	18.00
<b>Total</b>	<b>58.79 gm/liter</b>
pH (at 25°C )	6.5

\*\*Formula adjusted standard to suit the performance parameter

#### Principle And Interpretation

*Xanthomonas campestris pv.campestris* and *Xanthomonas campestris pv. armoraciae*, the causal agent of black rot of crucifers, is a bacterial pathogen of significant economic importance. The disease occurs in many countries and regions. It produces a range of extracellular enzymes (including proteases, pectinases and cellulases) and extra cellular polysaccharide (EPS), which collectively play essential roles in pathogenesis (2).

The original medium was developed by Chang (1) for the isolation of *Xanthomonas campestris* which was then modified by lowering the pH of the medium using additional potassium dihydrogen phosphate.

The medium contains soluble starch which helps in the detection of starch hydrolyzing organisms. Starch is a reserve sugar for the plant, hence phyto pathogen utilizes starch from plant as a energy source. Bacteria which show a clear zone around the growth produce the exoenzyme amylase which cleaves the starch into di- and monosaccharides.

Soya peptone and tryptone provide nitrogenous compounds, carbon, vitamin B complex and trace ingredients. The inorganic phosphates in the medium serve as buffers. Amino acids provide necessary growth factor for bacterial population. Glucose acts as a alternate carbon source for the organism.

#### Directions

Suspend 58.79 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121) for 15 minutes. Cool to 45-50C and aseptically add the rehydrated content of one vial of NNB supplement (PHS1015). Mix well and pour into Petri plate .

Note : Store plates for 4 days at 4C to improve visibility of starch hydrolysis.

### Quality Control

**Appearance:**

Cream to yellow coloured, homogeneous, free flowing powder.

**Gelling:**

Firm, comparable with 1.8% Agar gel.

**Colour and Clarity of prepared medium:**

Yellow coloured, opalescent gel with white precipitate forms in Petri plates

**Reaction:**

Reaction of 5.88 % w/v aqueous solution is pH 6.5 ± 0.2 at 25°C.

**Cultural Response:**

Cultural characteristics observed after an incubation at 25-30°C for 5-7 days.

Organism (ATCC)	Growth	Colony Characteristics
<i>Xanthomonas campestris</i> pv <i>campestris</i>	luxuriant	Yellow, Mucoïd colonies surrounded by a zone of starch hydrolysis
<i>Xanthomonas campestris</i> pv. <i>armoraciae</i>	luxuriant	Yellow, Mucoïd colonies surrounded by a zone of starch hydrolysis
<i>Escherichia coli</i>	inhibited	-
<i>Staphylococcus aureus</i>	inhibited	-
<i>Saccharomyces cerevisiae</i>	inhibited	-

### Storage and Shelf Life

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

### Further Reading

- Chang, C.J., Donaldson, R., Crowley, M, and Pinnow, D. 1991. A new semiselective medium for the isolation of *Xanthomonas campestris* pv *campestris*. *Phytopathology* 81: 449453.
- Dow, J. M., and Daniels, M. J. 1994. Pathogenicity determinants and global regulation of pathogenicity of *Xanthomonas campestris* pv. *campestris*. Pages 29-41 in: *Current Topics in Microbiology and Immunology*, Vol. 192: *Bacterial Pathogenesis of Plants and Animals*. J. L. Dangl, ed. Springer-Verlag, Berlin.

### Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performs parameters.