

Technical Information

Phyto Buffered Starch Agar Base

Product Code: PHM1002

Application: Selective medium for the detection of *Xanthomonas campestris pv.campestris* and *Xanthomonas campestris pv. armoraciae* in brassica.

Composition**

Ingredients	Grams/Litre
Yeast extract	0.10
Soluble starch	25.00
Di-potassium hydrogen phosphate	0.80
Potassium nitrate	0.50
Magnesium sulphate anhydrous	0.488
Methyl green	0.015
Agar	15.00
Final pH (at 25°C)	6.8

**Formula adusted standard to suit the performance parameter

Principle And Interpretation

This medium was used to identify the strain causing black rot in brassica. This medium is formulated as developed by D. Fieldhouse and M. Sasser. This medium is sensitive medium for the growth of *Xanthomonas campestris pv. campestris* and *Xanthomonas campestris pv. armoraciae*, the appearance of this colonies are greenish blue, translucent, mucoid and gives starch hydrolysis.

Yeast extract serves as a nutritional factor. Phosphates buffer the medium. Potassium nitrate serves a nitrogen source. Starch hydrolysis is detected by a visible zone of hydrolysis. The antibiotics in the supplement serves a selective agent.

Directions

Suspend 42.21 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and aseptically add the rehydrated contents of one vial of NMPT supplement (PHS002). Mix well and pour into Petri plates.

Note: Store plates for 4 days at 4°C to improve visibility of starch hydrolysis.

Quality Control

Appearance

Light yellow to light green coloured, homogeneous, free flowing powder.

Gelling

Firm, comparable with 1.5% Agar gel.

Colour and Clarity of prepared medium

Yellow coloured, with greenish tinge opalescent gel forms in Petri plates.

Cultural Response

Cultural characteristics observed with added NMPT Supplement, after an incubation at 30-32°C for 5-6 days.

Organism (ATCC)	Growth	Colony Characteristics	Starch Hydrolysis
<i>Xanthomonas campestris</i> <i>pv. Campestris</i>	luxuriant	pale green mucoid	+
<i>Xanthomonas campestris</i> <i>pv. Armoraciae</i>	luxuriant	pale green mucoid	+
<i>Saccharomyces cerevisiae</i> (9763)	inhibited	-	-

Key: Small zone surrounded the colony

Storage and Shelf Life

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Further Reading

- Yuen, G.Y., Alvarez, A.M., Benedict, A.A., and Trotter, K.J. 1987. Use of monoclonal antibodies to monitor the dissemination of *Xanthomonas campestris pv. campestris*. *Phytopathology* 77:366-370

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.