

Technical Information

Phyto Casein Hydrolysate Yeast Extract Agar Base

Product Code: PHM1004

Application: Semi selective medium for the detection of *Clavibacter michiganensis subsp michiganensis* on seeds of tomato.

Composition**

| Ingredients | Grams/Litre |
|-------------------------------|-------------|
| Yeast extract | 2.00 |
| Casein hydrolysate | 4.00 |
| Glucose, anhydrous | 10.00 |
| Magnesium sulphate, anhydrous | 0.15 |
| Boric acid | 1.00 |
| Ammonium chloride | 1.00 |
| Tris HCl | 1.20 |
| Agar | 18.00 |
| Final pH (at 25°C) | 7.4 |

**Formula adjusted standard to suit the performance parameter.

Principle And Interpretation

Clavibacter michiganensis subsp. michiganensis is a Gram-positive bacterium that causes wilting and cankers, leading to severe economic losses in commercial tomato production worldwide. Bacterial canker is one of the most destructive tomato diseases. Symptoms include necroses on leaves, stems and fruits during flowering of the plants.(1) It can be recovered on any general microbiological media such as nutrient agar, however specialized medium such as Phyto Casein Hydrolysate Yeast Extract Agar Base are more advantageous.(2). It is comparatively a less selective medium.

Casein hydrolysate and yeast extract supplies nitrogenous compounds. Glucose serves as a source of carbohydrate. Boric acid helps in selectivity of the medium. Tris HCL helps in buffering of the medium.

Directions

Suspend 37.35 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize both the solutions separately by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50 °C. Aseptically add the contents of one vial of CNP supplement (PHS003). Mix well and pour into sterile Petri plates

Quality Control

Appearance

Cream to yellow coloured, homogeneous, free flowing powder.

Gelling:

Firm, comparable with 1.8% Agar gel.

Colour and Clarity of prepared medium:

Yellow coloured opalescent gel with white precipitate forms in Petri plates

Reaction:

Reaction of 3.77% w/v aqueous solution is pH 7.4 at 25°C

Cultural Response:

Cultural characteristics observed with added CNP supplement (PHS1003), after an incubation at 25-30°C for 6-7 days.

| Organism (ATCC) | Growth | Colony characteristics |
|----------------------------------|-----------|----------------------------|
| <i>Clavibacter michiganensis</i> | Luxuriant | Yellow, mucoid, glistening |
| <i>subsp michiganensis</i> | | |
| <i>Escherichia coli</i> (25922) | inhibited | |

Storage and Shelf Life

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Further Reading

- Bradbury, J.F.(1986) Guide to plant pathogenic bacteria. CAB International Wallingford, UK.
- Fatmi, M., Schaad, N.W. 1988. Semiselective agar medium for isolation of *Clavibacter michiganensis subsp. michiganensis* from naturally infected tomato seeds.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performs parameters.