

Technical Information

Gresshoff & Doy (DBM2) Macroelements Solution (10X)

Product Code: PL1005

Composition**

Ingredients mg/Litre

Potassium nitrate	1000.00
Ammonium nitrate	1000.00
Calcium nitrate	241.11
Magnesium sulphate	17.09
Potassium phosphate monobasic	300.00
Potassium chloride	65.00
TOTAL gm/100ml	2.62

Directions

Measure out approximately 600ml of distilled water. While stirring add 100ml of macroelements stock solution. Rinse the original bottle with a small volume of distilled water to remove traces of the solution. Add desired heat stable supplements prior to autoclaving. Adjust the medium to the desired pH using 1N HCl/NaOH. Make up the final volume to 1000ml with distilled water. Sterilize the medium by autoclaving at 15 lbs or 121°C for 15 minutes. Cool the autoclaved medium to 45°C before adding the filter sterilized heat labile supplements. Dispense the desired amount of medium in sterile culture vessels

Principle and Interpretation

Gresshoff & Doy (DBM2) macroelements (10X) is a filter sterilized stock solution. The stock contains macroelements as described by Gresshoff & Doy. Add 100ml per litre of macroelements stock solution to the prepared medium to achieve the desired concentration.

Quality Control

Appearance

Colourless to yellow, clear solution.

pH

4.2 ±0.5 of 100ml macroelements stock solution.

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response :

Cultural condition :

Incubation period	: 5 weeks
Relative humidity	: 60% ± 2%
Temperature	: 22°C ± 2°C
Photoperiod (D: N) in hours	: 16:8

Cell Line	Type of Culture	Results
<i>Musa</i> species	Shoot culture	No structural deformity observed No necrotic tissues, Actively growing shoots, No toxicity to shoots
<i>Daucus</i> species	Callus culture	No necrotic tissues, Actively growing callus, No toxicity to callus

[The medium is prepared as per direction. The growth promoting activity of this macroelements stock solution is evaluated using two plant species viz. *Musa* species and *Daucus* species through three passages. Plant growth hormones (e.g. 2,4-D, NAA, Kinetin and 6-BAP) are added in suitable combinations and concentrations.]

Further Reading

1. Gresshoff P.M. & Doy C.H., Z. Pflanzenphysiol., (1974), 73, 132 - 141

Storage Temperature and Shelf Life

Store the macroelements stock solution at 2-8°C away from direct light. If possible, the entire content of each bottle should be used immediately, aseptically after opening or else the unused portion should be stored at 2-8°C. Avoid contaminating the liquid. Use before the expiry date.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.