



# **Product Specification**

cdhfinechemical.com

## **Technical Information**

## Earle's Balanced Salt Solution 1X

With Sodium bicarbonate and Phenol red Without Calcium and Magnesium

#### Product Code: TL2110

Application: All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation. The function of a salt solution is:

• To maintain the medium within physiological pH range.

• To maintain intracellular and extra cellular osmotic balance.

Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism.

Earle's balanced salt solution is designed to equilibrate with a 5% CO2 in air mixture. TL2110 is Earle's

balanced salt solution with sodium bicarbonate and phenol red hence requires the cells to be grown in a 5% CO<sub>2</sub> environment. It does not contain calcium and magnesium.

Composition**	
Ingredients	mg / Litre
INORGANIC SALTS	
Potassium chloride	0.400
Sodium bicarbonate	2.200
Sodium chloride	6.800
Sodium dihydrogen orthophosphate	0.122
OTHERS	
Dextrose anhydrous	1.000
Phenol red sodium salt	0.011

### **Quality control**

#### Appearance

Red colored, clear solution **pH** 7.00 - 7.60 **Osmolality in mOsm/Kg H2O** 265.00 -305.00 **Sterility** No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification. **Toxicity test** Passes **Endotoxin Content** NMT 1EU/ml





# **Product Specification**

cdhfinechemical.com

### Storage and Shelf Life

Store at 15-30 °C away from bright light. Shelf life is 24 months. Use before expiry date given on the product label.

#### **Disclaimer**:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.