



# **Product Specification**

cdhfinechemical.com

## **Technical Information**

## Hanks' Balanced Salt Solution 1X

With Phenol red

### Without Calcium, Magnesium sulphate and Sodium bicarbonate

## Product Code: TL2108

**Application:** All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation. The function of a salt solution is:

•To maintain the medium within physiological pH range.

To maintain intracellular and extra cellular osmotic balance.

• Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism. Hanks' balanced salt solution is designed to equilibrate with air, hence does not require CO<sub>2</sub> air mixture.

TL2108 is Hanks' balanced salt solution with phenol red. It is designed for use with cells maintained in less CO<sub>2</sub> environment or CO<sub>2</sub> free environment. It does not contain calcium, magnesium sulphate and sodium bicarbonate.

Composition**	
Ingredients	mg / Litre
INORGANIC SALTS	
Disodium hydrogen phosphate anhydrous	48.000
Potassium chloride	400.000
Potassium dihydrogen orthophosphate	60.000
Sodium chloride	8000.000
OTHERS	
D-Glucose	1000.000
Phenol red sodium salt	11.000

### Methodology

1. Add 4.7ml of 7.5% sodium bicarbonate solution (TCL1013) in 1X balanced salt solution prior to use.

## Material requied but not provided

7.5% sodium bicarbonate solution (TCL1013)





# **Product Specification**

cdhfinechemical.com

## Quality control

#### Appearance

Colourless, clear solution

#### рΗ

6.3-6.90

#### Osmolality in mOsm/Kg H<sub>2</sub>O

250.00 -290.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Toxicity test

Passes

#### Endotoxin Content

NMT 1EU/ml

## Storage and Shelf Life

Store at 15-30 °C away from bright light. Shelf life is 24 months. Use before expiry date given on the product label.

#### **Disclaimer**:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.