



# **Product Specification**

cdhfinechemical.com

## **Technical Information**

### MicSera Minimum Essential Medium Eagle

w/ Earle's Salts, NEAA and Sodium bicarbonate w/o L- Glutamine

1X Liquid Cell Culture Medium requiring reduced serum Supplementation

Product Code: RSL1008

Application:- MicSera media are based on the classical formulations supplemented with insulin, transferrin and other advanced nutrients. The additional nutrients help in reducing the percentage of serum required to grow most of the common cell lines. The percentage of serum reduction may vary with type of cell line used. For nonfastidious cell lines serum can be reduced from 10% to as low as 1%. For fastidious cell lines serum usage can be reduced from 10% to 2.5%. MicSera medium can be used without prior adaptation and sub cultured using normal procedures. Reduced serum supplementation improves the reproducibility of experimental results by decreasing the variability caused due to undefined serum constituents. It also facilitates down regulation process in bioassays and in purification process of culture products. MicSera Minimum Essential Medium Eagle is medium with Earle's salts, NEAA and Sodium bicarbonate. It does not contain L-glutamine. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

Composition**	
Ingredients	mg/Litre
INORGANIC SALTS	
Calcium chloride dihydrate	265.000
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium bicarbonate	2200.000
Sodium chloride	6800.000
Sodium dihydrogen phosphate anhydrous	122.000
AMINO ACIDS	
Glycine	7.500
L-Alanine	8.900
L-Arginine hydrochloride	126.000
L-Asparagine monohydrate	15.000
L-Aspartic acid	13.300
L-Cystine dihydrochloride	31.300
L-Glutamic acid	14.700
L-Histidine hydrochloride monohydrate	42.000
L-Isoleucine	52.000
L-Leucine	52.000
L-Lysine hydrochloride	72.500
L-Methionine	15.000
L-Phenylalanine	32.000
L-Proline	11.500
L-Serine	10.500
L-Threonine	48.000
L-Tryptophan	10.000
L-Tyrosine disodium salt	51.900
L-Valine	46.000





# **Product Specification**

cdhfinechemical.com

VITAMINS	VI	T	AΓ	И	ΙN	S
----------	----	---	----	---	----	---

1.000 Choline chloride D-Ca-Pantothenate 1.000 Folic acid 1.000 Nicotinamide 1.000 Pyridoxal hydrochloride 1.000 Riboflavin 0.100 Thiamine hydrochloride 1.000 i-Inositol 2.000

#### OTHERS

D-Glucose 1000.000
Growth Supplement mix Proprietary
Phenol red sodium salt 11.000

### Methodology

1. Add 20ml of 200mM L-glutamine (TCL1012) or MiGlutaXL supplement (TCL1030) for 1 litre of medium.

#### Recommendations for use with MicSera Media:

- 1. MicSera media have been optimized at 2.5% serum concentration for a broad range of cell culture applications. Recommended concentrations of serum using MicSera media ranges from 1-5%. However the concentration of serum used may need to be adjusted for specific cell types or applications to achieve optimal results. Titration of FBS concentration is recommended to determine maximum serum reduction.
- 2. MicSera media are provided as 1X solutions and need to be supplemented with 4mM Glutamine and required amount of reduced serum.
- 3. In case of antibiotics being used to control contamination, it is recommended to reduce the amount of antibiotics in proportion to the amount of serum reduced.

## Material required but not provided

L-Glutamine solution 200mM (TCL1012) MiGlutaXL Supplement (TCL1030)

Fetal Bovine Serum (BA2112/BA12432)

## **Quality control**

#### Appearance

Orangish red colored, clear solution.

### рΗ

7.00 -7.60

#### Osmolality in mOsm/Kg H2O

280.00 -320.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

#### Cultural Response

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control medium through minimum three subcultures.

#### Endotoxin Content

NMT 5EU/ml

## Storage and Shelf Life

- Store at 2-8°C away from bright light.
- Shelf life is 12 months.
- Use before expiry date given on the product label.





# **Product Specification**

cdhfinechemical.com

### Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.