

Technical Information

Heart Infusion Agar, Modified

Product Code: DM 1073F

Application: - Heart Infusion Agar, Modified is recommended in the isolation and cultivation of diarrheagenic *Escherichia coli* in accordance with FDA BAM, 1998.

Composition**

Ingredients	Gms / Litre
Heart muscle, infusion from	7.500
Soya peptone	10.000
Sodium Chloride	5.000
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

A liquid medium containing an infusion of meat was one of the first media used for the cultivation of bacteria. These infusion media need not be further supplemented by the addition of supplements incase of fastidious bacteria (1). Heart Infusion Agar, Modified is used in the isolation and cultivation of diarrheagenic *Escherichia coli* in accordance with FDA BAM, 1998 (2). *E.coli* is a Gram negative, facultatively anaerobic bacterium that is found as commensals in human intestine.

Heart Infusion Agar, Modified is cited in BAM for primary screening in the conventional biochemical screening and identification of diarrheagenic *Escherichia coli*. Enrichment of the sample on BHI and TP broths is recommended as the first step in primary screening of diarrheagenic *E. coli*. This may induce the growth and proliferation of other members of *Enterobacteriaceae* including non lactose fermenting strains of *E.coli*. Hence additional tests may need to be performed for isolation. Transfer suspicious colonies to TSI agar, Heart Infusion Agar, Modified slants (DM 1073F), Tryptone Broth, Arabinose Broth, and Urea Broth and incubate for 20 hr at 35°C. Organisms isolated on primary screening are processed for secondary screening and confirmed using genotypic, biochemical and serological reactions. On supplementation of blood, Heart Infusion Agar, Modified can be used to study haemolytic reactions (3). Heart Muscle, infusion from and soya peptone provide nutritional requirements for the pathogenic bacteria. Sodium chloride maintains the osmotic equilibrium of the medium.

Methodology

Suspend 37.5 grams of dehydrated media powder in 1000 ml distilled water. Mix thoroughly & heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 20 minutes. Cool to 50°C. Mix well and pour into sterile Petri plates.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity

Basal medium: Light amber coloured clear to slightly opalescent gel. After addition of 5% v/v sterile defibrinated blood: Cherry red coloured opaque gel forms in Petri plates.

Reaction

Reaction of 4.0% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH Range

7.10-7.50

Cultural Response

DM 1073: Cultural characteristics observed with added 5% w/v sterile defibrinated blood, after an incubation at 35-37°C for 18-48 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth w/o blood	Recovery w/o blood	Growth with blood	Recovery with blood	Haemolysis
Cultural Response						
<i>Neisseria meningitidis</i> ATCC13090	50-100	fair	40-50%	luxuriant	≥70%	none
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good	50-70%	luxuriant	≥70%	beta
<i>Staphylococcus epidermidis</i> ATCC 12228	50-100	good	50-70%	luxuriant	≥70%	none
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	fair-good	40-50%	luxuriant	≥70%	alpha
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	fair-good	40-50%	luxuriant	≥70%	beta

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and the prepared medium at 2 -8°C. Use before expiry date on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Hansen, N. H. 1962. J. Appl. Bacteriol., 25.
2. FDA, U.S. 1998. Bacteriological Analytical Manual. 8 ed. Gaithersburg, MD: AOAC International.
3. Downes, F.P. and Ito, K. 2001. Methods For The Microbiological Examination of Foods. APHA, Food 4 ed. Washington, D.C.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specification for identity and performance parameters.