

Technical Information

LPM Agar Base

Product Code: DM 2228

Application: - LPM (Lithium Phenylethanol Moxalactam) Agar is recommended for isolation and cultivation of *Listeria monocytogenes* from food and dairy products.

Composition**					
Ingredients	Gms / Litre				
Casein enzymic hydrolysate	5.000				
Peptic digest of animal tissue	5.000				
Beef extract	3.000				
Glycine anhydride	10.000				
Lithium chloride	5.000				
Sodium chloride	5.000				
Phenylethyl alcohol	2.500				
Agar	15.000				
Final pH (at 25°C)	7.3±0.2				
**Formula adjusted, standardized to suit performa	nce parameters				

Principle & Interpretation

L. monocytogenes is a gram-positive foodborne human pathogen responsible for serious infections in pregnant women that may ultimately result in abortion, stillbirth, birth of child with meningitis or primary bacteremia in adults and juveniles. Lee and McClain (1) developed LPM Agar, which is a modification of McBride Listeria Agar. It enhances the recovery of low numbers of L. monocytogenes from mixed microflora in samples. APHA also recommends this medium for food and dairy sample testing (2, 3).

The first reported foodborne outbreak of listeriosis was in 1985 (6) and since then, microbiological and epidemiological evidence from both sporadic and epidemic cases of listeriosis has shown that the principal route of transmission is via the consumption of foodstuffs contaminated with *Listeria monocytogenes* (7). The organism has been isolated from commercial dairy and other food processing plants, and is ubiquitous in nature, being present in a wide range of unprocessed foods and in soil, sewage and silage and river water (8). *Listeria* species grow over a pH range of 4.4 to 9.6, and survive in food products with pH levels outside these parameters (2). Motility is most pronounced at 20°C.

Peptic digest of animal tissue, casein enzymic hydrolysate and beef extract supply nitrogen sources, vitamins and minerals for the test organisms. Sodium chloride maintains the osmotic balance of the medium. Glycine anhydride improves recovery of Listeria. Lithium chloride, moxalactam and phenyl ethanol aids in suppression of both gram-positive and gram-negative organisms including *Staphylococcus*, *Proteus* and *Pseudomonas* species. *Listeria monocytogenes* show blue-green iridescence when examined with oblique transmitted light (4, 5). For liquid samples, pipette 25 ml sample in to 225 ml Listeria Enrichment Broth in 500 ml flask and mix well by shaking. For solid samples, weigh 25 grams into a stomacher bag and add 225 ml of enrichment broth. Mix the sample by stomaching for 2 minutes and incubate the mixture in the plastic bag at 30°C. After incubation for 24 hours and again at 48 hours, the enriched culture is streaked onto LPM Agar (DM2228). Incubate at 35°C for 48 hours and presumptive *Listeria* colonies are selected under 45°C transillumination.





Methodology

Suspend 50.5 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat to boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 12 minutes. Cool to 50°C Aseptically add rehydrated contents of 1 vial of Moxalactam Supplement (MS 2151). Mix well before pouring into sterile Petri plates.

Warning: Lithium chloride is harmful. Avoid bodily contact and inhalation of vapours. On contact with skin, wash with plenty of water immediately.

Quality Control

Appearance

Cream to light yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity

Light yellow coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 5.05% w/v aqueous solution at 25°C. pH: 7.3±0.2

pH Range

7.10-7.50

Cultural Response

DM 2228: Cultural characteristics observed with added Moxalactum Supplement (MS 2151) after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
Cultural Response Escherichia coli ATCC 25922	>=10³	inhibited	0%
Listeria monocytogenes ATCC 19111	50-100	good-luxuriant	>=50%
Listeria monocytogenes ATCC 19112	50-100	good-luxuriant	>=50%
Listeria monocytogenes ATCC 19117	50-100	good-luxuriant	>=50%
Pseudomonas aeruginosa ATCC 27853	>=10³	inhibited	0%
Staphylococcus aureus ATCC 25923	>=103	inhibited	0%

Storage and Shelf Life

Dried Media: Store below 10-30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label. **Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

Further Reading

- 1. Lee W. H. and McClain D., 1986, Appl. Environ. Microbiol., 52:1215.
- 2. Downes F. P and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
- 3. Marshall R. T., (Ed.), 1992, Standard Methods for the Examination of Dairy Products, 16th Ed., APHA, Washington, D.C.
- 4. Bearns R. E. and Girard K. F., 1959, Am. J. Med. Technol., 25:120.
- 5. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Yolken R. H., (Ed.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- 6. Wehr H. M., 1987, J. Assoc. Off. Anal. Chem., 70:769.
- 7. Bremer P. and Osborne C., 1995, J. Food Prot., 58:604.
- 8. Patel J. R., Hwang C. A, Beuchat L. R, Doyle M. P. and Brackett R. E., 1995, J. Food Prot., 58:244





Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specification for identity and performance parameters.

