

Technical Information

Listeria Selective Agar (Twin Pack)

Product Code: DM 1567

Application: - Listeria Selective Agar is recommended for cultivation and selective isolation of *Listeria* species from clinical specimens.

Composition**

Ingredients	Gms / Litre
Part A	-
Casein enzymic hydrolysate	10.000
Peptic digest of animal tissue	10.000
Dextrose	1.000
Sodium chloride	5.000
Thiaminium dichloride	0.005
Acriflavin hydrochloride (Trypaflavin)	0.010
Nalidixic acid	0.040
Agar	13.000
Part B	-
Potassium thiocyanate	37.500
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Listeria Selective Agar was proposed by Feindt (1) for the cultivation of *Listeria* species from clinical and non-clinical specimens. Obiger and Schonberg (2) reported the superiority of these media to yield *Listeria* from mix-infected specimens.

Casein enzymic hydrolysate, peptic digest of animal tissue supplies essential nutrients. Thiaminium dichloride is the vitamin B source added to improve the growth of *Listeria*. Thiocyanate and Nalidixic acid prevents gram-negative bacteria (3, 4). Bockemühl (5) reported suppression of Enterococci by combination of selective agents and acridine dyes. The combination of Acriflavin hydrochloride and Nalidixic acid was used by Ralovich et al (6) and Kampelmacher and Van Noorle Jansen (7) for the isolation of *Listeria*. Listeria Enrichment Broth can be further improved by adding Colimycin alongwith Nalidixic acid (8). The mix infected specimen is added directly to Listeria Enrichment Broth or subjected to cold enrichment and then cultured on Listeria Selective Agar. Haemolytic forms can be identified by inoculating Blood Agar (DM 1073).

Methodology

Suspend 39 grams of dehydrated powder media **Part A** and 37.5 grams of **Part B** in 1000 ml distilled water. Mix thoroughly & heat to boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Appearance

Part A: Cream to yellow homogeneous free flowing powder

Part B: White to cream homogeneous free flowing powder

Gelling

Firm, comparable with 1.3% Agar gel.

Reaction

Reaction of medium (3.9% w/v **Part A** + 3.75% w/v **Part B**) at 25°C. pH : 7.4±0.2

pH Range

7.20-7.60

Cultural Response

DM 1567: Cultural characteristics observed in presence of 10% Carbon dioxide (CO₂) after an incubation at 35-37°C for 48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Enterococcus faecalis</i> ATCC 29212	50-100	none-poor	<=10%
<i>Escherichia coli</i> ATCC 25922	>=10 ³	inhibited	0%
<i>Listeria innocua</i> ATCC 33090	50-100	luxuriant	>=50%
<i>Listeria ivanovii</i> ATCC 19119	50-100	luxuriant	>=50%
<i>Listeria monocytogenes</i> ATCC 19112	50-100	luxuriant	>=50%
<i>Listeria monocytogenes</i> ATCC 19118	50-100	luxuriant	>=50%

Storage and Shelf Life

Dried Media: Store dehydrated and prepared medium at 2-8°C. Use before expiry period on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Feindt E., 1972, Inuug. Diss., Würzburg.
2. Obiger G. and Schonberg A., 1973, Fleischwirtschaft, 10:1450.
3. Lebnert C., 1964, Arch. Exp. Vet. Med., 8: 891 and 1247.
4. B eerens H. and Tahon-Castel M.M., 1966, Ann. Inst. Pasteur, 111:90.
5. Bockemühl J., Seeliger H.P.R. and Kathke R., 1971, J. Med. Microbiol. Imm. 157:84.
6. Ralorich B., et al, 1971, Zbl. Bakt. I.Orig., 216:88.
7. Kampelmacher E.H. and Van Noorle-Jansen L.M., 1972, Zbl. Bakt. J.Orig., 221:139.
8. Grey M.L. et al, 1948, J. Bact., 55:471.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specification for identity and performance parameters.