

Bases / Media Supplements

# **Technical Information**

## **Tryptone Soya Yeast Extract Broth**

# Product Code: DM 2263

Application: - Recommended for confirmation of *Listeria* in Henry's light.

Composition**		
Ingredients	Gms / Litre	
Tryptone	17.000	
Soya peptone	3.000	
Sodium chloride	5.000	
Dipotassium hydrogen phosphate	2.500	
Dextrose (Glucose)	2.500	
Yeast extract	6.000	
Final pH ( at 25°C)	7.3 <u>+</u> 0.2	
**Formula adjusted, standardized to suit performance	e parameters	

### **Principle & Interpretation**

Tryptone Soya Yeast Extract Broth is formulated as per APHA (1) for the isolation and cultivation of *Listeria monocytogenes* from foods. ISO Committee (2) has recommended for the cultivation and maintenance of a wide variety of heterotrophic microorganisms (3).

Tryptone and soya peptone provide nitrogeneous and carbonaceous compounds, long chain amino acids and other essential nutrients. Dextrose is the energy source. Dipotassium hydrogen phosphate acts as buffering system to control pH. Yeast extract is the rich source of vitamin B complex.

According to FDAs enrichment procedure (4) for isolation of *Listeria monocytogenes* from dairy products, the sample to be tested is inoculated in enrichment broth and incubated at 30°C for 24-48 hour

### Methodology

Suspend 36 grams in 1000 ml purified / distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

**Note:** Just before use, heat the test tubes on a boiling water-bath or in flowing steam for 10 min, and then cool rapidly to about 45°C.

# Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder.

#### Colour and Clarity

Yellow coloured clear solution in tubes.

#### Reaction

Reaction of 3.6% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH Range

7.10-7.50

#### Cultural Response

DM2263: Cultural characteristics observed after an incubation at 30-37°C for 24-48 hours.





Bases / Media Supplements

Organism	Inoculum	Growth	
	(CFU)		
Listeria monocytogenes ATCC 19117	50-100	good-luxuriant	
Listeria monocytogenes ATCC 19111 (00020*)	50-100	good-luxuriant	
Listeria monocytogenes ATCC 19118	50-100	good-luxuriant	

Key : \*Corresponding WDCM numbers.

### Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 15-25°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle inorder to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

### Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

### **Disclaimer :**

• User must ensure suitability of the product(s) in their application prior to use.

• The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate

- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.

• Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performances parameters.

