

Technical Information

Trptone Agar

Product Code: DM 2365

Application: - Trptone Agar is used as a general purpose medium for growth of non-fastidious microorganisms.

Composition**

Ingredients	Gms / Litre	
Casein enzymic hydrolysate	10.000	
Sodium chloride	8.000	
Agar	15.000	
**Formula adjusted, standardized to suit performance		

Principle & Interpretation

Tryptone Agar is a general purpose nutritious medium for growth of non-fastidious microorganisms. Tryptone Agar was developed by Vera (1) for the accurate differentiation and identification of aerobes and anaerobes by means of motility and fermentation reactions. It is recommended for Clostridia, *Bacillus* species, Micrococci, enteric bacilli and other non-fastidious organisms (2). This is also an excellent medium for the maintenance for both aerobic and anaerobic cultures. Viability in this medium is greater than in any other broth medium or slant culture.

Casein enzymic hydrolysate provides essential growth nutrients to support the growth of organisms. Sodium chloride buffers the medium.

Methodology

Suspend 33.0 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured clear to slightly opalescent gel forms in Petri plates.

Cultural Response

DM 2365 Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
Enterobacter aerogenes ATCC 13048	50-100	luxuriant	>=70%
Escherichia coli ATCC 25922	50-100	luxuriant	>=70%
Pseudomonas aeruginosa ATCC 27853	50-100	luxuriant	>=70%
Salmonella Enteritidis ATCC 13076	50-100	luxuriant	>=70%
Staphylococcus aureus ATCC 25923	50-100	luxuriant	>=70%





Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Further Reading

- Vera, 1944, J. Bacteriol., 47:455.
- 2. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

