

## Technical Information

### MUG Lauryl Sulphate Broth, Modified

#### Product Code: DM 2046I

**Application:** - MUG Lauryl Sulphate Broth, Modified is recommended as a selective enrichment medium for enumeration of presumptive *Escherichia coli* and other coliforms from milk and milk products.

#### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	20.000
Lactose	5.000
Dipotassium hydrogen phosphate	2.750
Monopotassium phosphate	2.750
Sodium chloride	5.000
4-Methylumbelliferyl β-D-glucuronide (MUG)	0.100
Sodium lauryl sulphate	0.100
Tryptophan	1.000
Final pH ( at 25°C)	6.8±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

MUG Lauryl Sulphate Broth, Modified is recommended by the ISO Committee (6) for enumeration of presumptive *E.coli* and other coliforms from milk and milk products. The coliform group consists of several genera of bacteria belonging to the family *Enterobacteriaceae*. Lauryl Sulphate Broth formulated by Mallmann and Darby (1) is recommended by APHA for the detection and enumeration of coliform organisms in foods, water and waste water (2, 3). MUG is added in Lauryl Sulphate Broth as the fluorogenic compound which permits the rapid detection of *Escherichia coli* when observed under UV light where further confirmation is not required (2,4).

*E.coli* possesses the enzyme beta-glucuronidase and is capable of cleaving the fluorogenic substrate 4-Methylumbelliferyl beta-D-glucuronide (MUG) with the corresponding release of the fluorogen.

MUG detects anaerogenic strains, which may not be detected in the conventional procedure. Feng and Hartman (5) used MUG- containing medium for studying beta-glucuronidase activity and found that *E.coli* has 96-100% activity, while *Salmonella* species has 17% and *Shigella* species 40% activity and other genera were negative. For weakly positive strains, incubation should be carried out overnight. Robison (4) achieved no false negative results and about 5% false positive results.

Casein enzymic hydrolysate provides nutrients while lactose serves as energy source. Sodium lauryl sulphate prevents many organisms other than coliforms. 4-methylumbelliferyl-beta-D-glucuronide is hydrolyzed by the enzyme beta-glucuronidase to yield 4-methylumbelliferone, a fluorescent end product.

Inoculate 10 ml of the test specimen into three tubes each of single strength and double strength. Incubate the tubes at 35°C for 24 hours. Observe for opacity and gas formation. For confirmation of presumptive *E.coli*, observe for fluorescence and perform indole reaction using Kovacs Reagent (6).

#### Methodology

Suspend 36.7 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat if necessary to dissolve the medium completely. Dispense into tubes with inverted Durham's tubes as required, taking into account the volume of the sample to be tested. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.



Dehydrated Culture Media  
Bases / Media Supplements

## Quality Control

### Appearance

Cream to yellow homogeneous free flowing powder

### Colour and Clarity

Light amber coloured clear solution without any precipitate

### Reaction

Reaction of 3.67% w/v aqueous solution at 25°C. pH : 6.8±0.2

### pH Range

6.60-7.00

### Cultural Response

DM 2046I: Cultural characteristics observed after an incubation at 30°C for 24 hours.

Organism	Inoculum (CFU)	Growth	Fluorescence under uv	Indole production
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	negative	negative reaction
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	positive	positive reaction, red ring at the interface of the medium

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Mallmann W. L. and Darby C. W., 1941, Am. J. Public Health, 31:127.
2. Downes F. P and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
3. Greenberg A. E., Trussell R. R. and Clesceri L. S., (Eds.), 1985, Standard Methods for the Examination of Water and Wastewater, 16th Ed., APHA, Washington, D.C.
4. Robison, 1984, Appl. Environ. Microbiol., 48:285.
5. Feng P. C. S. and Hartman P. A., 1982, Appl. Environ. Microbiol., 43 : 1320.
6. International Organization for Standardization (ISO), ISO 11866-2:1997 (E)

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

