

Technical Information

M-PA Agar

Product Code: DM 2121A

Application: - M-PA Agar is recommended for the selective detection and isolation of *Pseudomonas aeruginosa* by membrane filter technique.

Composition**

Ingredients	Gms / Litre
Yeast extract	2.000
L-Lysine hydrochloride	5.000
Sodium chloride	5.000
Xylose	1.250
Sucrose	1.250
Lactose	1.250
Ferric ammonium citrate	0.800
Sodium thiosulphate	5.000
Phenol red	0.080
Magnesium sulphate	1.500
Kanamycin	0.008
Nalidixic acid	0.037
Agar	12.000
Final pH (at 25°C)	7.2±0.1

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

The MPN (Most Probable Number) technique results in satisfactory recovery levels of *Pseudomonas aeruginosa*, but is not usable for the testing of large-volumes water samples and they also lack precision for the recovery of *Pseudomonas aeruginosa*.

Levin and Cabelli devised M-PA Agar as a selective membrane filter medium for *Pseudomonas aeruginosa* (1). Many of the filter media used for the recovery of *P.aeruginosa* lack specificity and are of limited value when a large heterogeneous microbial flora is present in the water samples. M-PA Agar incorporates four antibiotics kanamycin, nalidixic acid, sulfapyridine and cycloheximide which render the medium moderately selective. The original medium was modified by raising the pH (2) and altering the content or concentration of ingredients (3). This media is included in part 914 C, Membrane Filter Technique for *P.aeruginosa* (Tentative) in the 16th / 19th Edition of Standard Methods for the Examination of Water and Waste water (4).

Brodsky and Ciebin (5) further modified these media by eliminating Sulfapyridine and Cycloheximide. This modified M-PA Agar helps in enumeration of *Pseudomonas aeruginosa* within 24 hours of incubation at 41.5°C compared to 72 hours for a presumptive M-PN test.

Yeast extract, lysine and carbohydrates provide nitrogenous compounds, energy sources and vitamins required for bacterial metabolism. Sodium chloride maintains osmotic equilibrium. Inorganic salts provide essential ions. Kanamycin inhibits protein synthesis in gram-positive organisms (6). Cycloheximide (MS1202) inhibits fungal flora. Nalidixic acid blocks replication of susceptible gram-negative bacteria (6). Phenol red is the pH indicator which turns yellow under acidic conditions due to fermentation of the carbohydrates.

After filtration of the water sample through a sterile 0.45 µm gridded filter, place the membrane filter on the surface of plates of M-PA Agar taking care to avoid the entrapment of bubbles between the agar and filter surface. Incubate for 24 hours at 41.5±0.5°C in an aerobic atmosphere. Optimal colony density on membrane filters is 20-200 colonies. All colonies on the filter Surface inoculate the above diluted cultures and incubate at 28°C for 48 hours to determine MIC value of fluconazole. are counted when the density is 2 or fewer per square; the average of 10 squares is determined when the count is 3-10 colonies per square and the average of 5 squares is determined when the count is 10-20 colonies per square. The average count per square is multiplied by 100 times the reciprocal of the dilution to give colonies per ml.

Methodology

Suspend 35.18 grams of dehydrated powder media in 1000 ml of distilled water. Mix thoroughly & heat with agitation and boil for 1 minute to dissolve the medium completely. DO NOT AUTOCLAVE. Cool to 45-50°C and pour into sterile Petri plates. Use the medium within 1 week of preparation.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.2% Agar gel.

Colour and Clarity

Orange red coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.52% w/v aqueous solution at 25°C. pH : 7.2±0.1

pH Range

7.10-7.30

Cultural Response

DM 2121A: Cultural characteristics observed after an incubation at 41.5 ± 0.5°C for up to 72 hours.

Organism	Inoculum (CFU)	Growth	Colour of medium
<i>Escherichia coli</i> ATCC 25922	>=10 ³	inhibited	-
<i>Klebsiella pneumoniae</i> ATCC 13883	>=10 ³	inhibited	-
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	good-luxuriant	pink
<i>Salmonella Typhi</i> ATCC 6539	>=10 ³	inhibited	-
<i>Staphylococcus aureus</i> ATCC 25923	>=10 ³	inhibited	-

Storage and Shelf Life

Dried Media: Store dehydrated powder and prepared medium at 2 - 8°C. Use before expiry date on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Levin M. A. and Cabelli V. J., 1972, Appl. Microbiol., 24:864.
2. Carson L. A., Peterson M, J., Favero M. S., Doto I. L., Collins D. E. and Levin M. A., 1975, Appl. Microbiol., 30:935.
3. Dutka B. J. and Kwan K. K., 1977, Appl. Environ. Microbiol., 33:240.



Dehydrated Culture Media
Bases / Media Supplements

4. Greenberg A. E. , Trussell R. R. and Clesceri L. S., (Eds.), 1985, Standard Methods for the Examination of Water and Wastewater, 16th / 19th Ed., APHA, Washington, DC.
5. Brodsky M. H. and Ciebin Y. B. W., 1978, Appl. Environ., Microbiol,36:36.
6. Estevez R. A., 1984, Bacteri

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specification for identity and performance parameters.

