

# **Technical Information**

## **Antibiotic Assay Medium No.1**

## Product Code: DM 1003U

Application: - Antibiotic Assay Medium No.1 is recommended for microbiological assay of ß-lactam and other antibiotics in accordance with United States Pharmacopoeia.

## Composition\*\*

•		
Ingredients	Gms / Litre	
Peptone	6.000	
Pancreatic digest of casein	4.000	
Yeast extract	3.000	
Beef extract	1.500	
Dextrose	1.000	
Agar	15.000	
pH after sterilization	6.6±0.1	
**Formula adjusted, standardized to suit perform	ance parameters	

### Principle & Interpretation

This medium is also recommended as inoculum and maintenance medium for different test organisms for antibiotic assays. Composition of this medium is in accordance with US Pharmocopoeia (1) and is recommended by FDA (2) and identified numerically with the name assigned by Grove and Randall (3).

Essential nutrients, vitamins, mineral, trace elements and growth factors are supplied by peptone, pancreatic digest of casein, yeast and beef extract. Dextrose in the medium serves as the carbon source for stimulating the growth of the test microorganism. Agar supplies excellent medium for antibiotic diffusion and gives well defined zones of inhibition.

Freshly prepared plates should be preferably used for assaying antibiotics. Test organisms is inoculated in sterile seed agar pre-cooled to 40-45°C and spread evenly over the surface of solidified base agar. All conditions in the microbiological assay must be controlled carefully. One of the critical and important step for obtaining good results is use of appropriate standard culture media.

## Methodology

Suspend 30.5 grams of dehydrated media in 1000 ml distilled water. Mix thoroughly & heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Advice: Recommended for the microbiological assay of Bacitracin, Nafcillin, Cephalothin,Cephaperin,Cloxacillin, Novobiocin, Penicillin-G

## Quality Control

#### Appearance

Cream to yellow coloured homogeneous free flowing powder

#### Gelling

Firm, comparable with 1.5% Agar gel

#### Colour and Clarity

Yellow coloured slightly opalescent gel forms in Petri plates.

### pH Range

6.50-6.70





#### **Cultural Response**

DM 1003U: Cultural characteristics observed after an incubation at specified temperature and period.

Organism	Inoculum (CFU)	Growth	Recovery	Inoculum medium	Incubation temperature / period	Antibiotics assayed
Bordetella bronchiseptica ATCC 4617	50-100	good-luxuriant	>=50%	Colistimethate sodium, Colistin, Polymyxin B	32-35°C/ 24 hours	
Escherichia coli ATCC 10536	50-100	luxuriant	>=70%	Chloramphenicol	32-35°C/ 24 hours	
Klebsiella pneumoniae ATCC 10031	50-100	good-luxuriant	>=50%	Capreomycin, Dihydrostreptomycin, Neomycin, Streptomycin, Troleandomycin	36-37.5°C/	
Micrococcus luteus ATCC 9341	50-100	luxuriant	>=70%	Erythromycin	32-35°C/ 24 hours	
Micrococcus luteus ATCC 10240	50-100	good-luxuriant	>=70%	Bacitracin	32-35°C/ 24 hours	Bacitracin
Pseudomonas aeruginosa ATCC 25619	50-100	luxuriant	>=70%	Carbenicillin	36-37.5°C/ 16-24 hours	
Staphylococcus epidermidis ATCC 12228	50-100	good-luxuriant	>=70%	Gentamicin, Netilmicin, Neomycin, Novobiocin, Paromomycin, Sisomycin	32-35°C/ 24 hours	Novobiocin
Staphylococcus aureus ATCC 29737	50-100	luxuriant	>=70%	Amikacin, Cephalothin, Cephaperin, Chlortetracycline, Cloxacillin, Cycloserine, Demeclocycline, Kanamycin, Methacycline, Nafcillin, Penicillin-G, Rolitetracycline, Tetracycline, Tobramycin, Tylosin	32-35°C/ 24 hours	Cephalothin, Cephaperin, Cloxacillin, Nafcillin, Penicillin-G,

# Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use freshly prepared medium. Use before expiry date on label. **Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.





# **Further Reading**

- 1. United States Pharmacopoeia USP 2011, US Pharmacopoeial Convention, Inc., Rockville, MD.
- 2. Tests and Methods of Assay of Antibiotics and Antibiotic containing Drugs, FDA, CFR, 1983 Title 21, Part 436, Subpart D, Washington, D.C.: U.S. Government Printing Office, paragraphs 436, 100-436, 106, p. 242-259, (April 1).
- 3. Grove and Randall, 1955, Assay Methods of Antibiotics Medical Encyclopaedia, Inc. New York

### Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performens parameters.

