

## Technical Information

### K Agar

**Product Code: DM 2752**

**Application:** - K Agar is recommended for isolation and cultivation of *Alicyclobacillus* in fruit juices in accordance Official Method of IFU.

### Composition\*\*

Ingredients	Gms / Litre
Yeast extract	2.500
Peptone	5.000
Glucose	1.000
Tween	1.000
Agar	15.000
Final pH ( at 25°C)	3.7±0.1

\*\*Formula adjusted, standardized to suit performance parameters

### Principle & Interpretation

Acidified environment are required to detect and isolate *A. acidoterrestris* (1, 2), therefore, K-Agar is recommended for detection of taint producing *Alicyclobacillus acidoterrestris* as per IFU (1) (standard IFU method No.12).

*Alicyclobacillus* are aerobic thermophilic non-pathogenic, spore forming bacteria that can survive the relatively mild pasteurization temperatures used for fruit juices and concentrates. Even very low numbers of *Alicyclobacillus* cause spoilage and off flavors in the beverages, damaging the brand (3). These bacteria are able to grow at pH values as low as 2.5 and also at elevated temperatures above 200°C. Their spores survive for long period in fruit concentrates and similar environments.

*A. acidoterrestris* is the most commonly occurring species that produce taints in juice and similar products, however other species may also produce taints.

Peptone and yeast extract acts as a source of nitrogen, amino acids, vitamins, and other essential growth requirements. Glucose acts as a carbon source. Polysorbate 80 serves as an additional source of growth factor and fatty acid. The low acidic pH (3.7) of medium obtained by addition of L-malic acid is inhibitory to several bacterial species.

K-Agar (when incubated at 45°C) supports the growth of predominant *A. acidoterrestris* and limited growth of other species. If the sample is filterable, filter it through 0.45µm membrane after pretreatment. Filter aseptically and transfer one membrane on to K- Agar and other membrane on YSG Agar (DM2753). Simultaneously, streak standard cultures of *Alicyclobacillus acidoterrestris* and *Alicyclobacillus acidocaldarius* for comparison. Incubate at 45±1°C for 2-5 days. Examine daily, for detailed procedure refer standard IFU method NO.12 (1)

### Methodology

Suspend 24.5 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to about 50°C and adjust pH 3.7 ± 0.1 with 25% L-malic acid.

### Quality Control

#### Appearance

Cream to yellow homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel.

### Colour and Clarity

Pale yellow coloured Clear to slightly opalescent gel forms in Petri plates.

### Reaction

Reaction of 2.45% w/v aqueous solution at 25°C. pH : 3.7±0.1

### pH Range

3.60-3.80

### Cultural Response

DM2752: Cultural characteristics observed after an incubation at 45- 46°C for 2 -5 days.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Alicyclobacillus acidocaldarius</i> ATCC 27009	50-100	Poor-good	30-40%
<i>Alicyclobacillus acidoterrestris</i> ATCC 49028	50-100	luxuriant	>=50%
<i>Alicyclobacillus acidocaldarius</i> ATCC 43030	50-100	Poor-good	30-40%
<i>Escherichia coli</i> ATCC 25922	>=10 <sup>3</sup>	inhibited	0%
<i>Staphylococcus aureus</i> ATCC 25923	>=10 <sup>3</sup>	inhibited	0%

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. IFU (2004). Standard IFU method No.12- Microbiological detection of Alicyclobacillus in fruit juices. International Federation of fruit juice producer, Paris
- 2 Matsubara et al. (2002). Alicyclobacillus acidiphilus sp. nov., a novel thermo-acidophilic-alicyclic fatty acid-containing bacterium isolated from acidic beverages. Int. J. Syst. Environment. Microbiol. 52, 1681-1685.
3. Baungart and Merve S. The Impact of Alicyclobacillus acidoterstris on the quality of Juices and Soft Drinks Fruit processing 7 : 251-254 (2000)

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.