

## Technical Information

### Modified McBride Listeria Agar Base

#### Product Code: DM 1891

**Application:** - Modified McBride Listeria Agar base is recommended for selective isolation and cultivation of *Listeria monocytogenes* from foodstuffs, clinical samples etc.

#### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	5.000
Peptic digest of animal tissue	5.000
Beef extract	3.000
Sodium chloride	5.000
Glycine anhydride	10.000
Lithium chloride	0.500
Phenyl ethanol	2.500
Agar	15.000
Final pH ( at 25°C)	7.3±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

The disease listeriosis is a frequent cause of abortions in cattle and sheep. In human, symptoms are manifested as septicemia, encephalitis and circulatory monocytosis (1). In women, *Listeria* may be harboured in the genital tract and may be transmitted transplacentally from infected amniotic fluid and vaginal discharge to infants (2). *Listeria* was first definitively described by Murray et al (3) in connection with an epizootic disease among laboratory-raised guinea pigs and rabbits. Since *Listeria* species are generally associated with animals, outbreaks of listeriosis are usually associated with contaminated dairy products. *Listeria* multiplies over a wide range of temperatures, from 3°C to 45°C, and over a pH range of 5.0 to 9.6. It also survives in food products with pH levels outside these parameters (4). Because of these properties, *Listeria* survives the various food processing techniques (5).

Modified McBride Listeria Agar Base differ from McBride Listeria Agar Base in the nutrient source available to *Listeria* species.

Casein enzymic hydrolysate and beef extract in the medium supply nitrogen, carbon, sulphur and trace nutrients required for the growth of *Listeria*. Phenyl ethyl alcohol is bacteriostatic for gram-negative bacteria as it selectively inhibits DNA synthesis (6). Sodium chloride helps to maintain the osmotic balance of the medium. Glycine inhibits certain gram-negative and gram- positive bacteria including *Escherichia coli* and *Enterococcus faecalis*, the common accompanying contaminants. Lithium chloride also has antibacterial activity. Further selectivity is achieved by the addition of McBride Listeria Supplement (MS2070). The detection of *L.monocytogenes* is greatly improved by pre-enrichment in liquid media either by one step or two steps. In one step method (7), infected material is inoculated directly in Listeria Selective Broth Base (DM1889), while in two steps method (8), infected material is inoculated in Listeria Enrichment Broth Base (UVM) (DM1890A) and incubated at refrigeration temperature of 4°C for few weeks (cold enrichment), as the organism has the ability to grow in low temperature. It is then inoculated in Fraser Secondary Enrichment Broth Base (DM2083), followed by plating onto selective agar such as Modified McBride Listeria Agar. The presumptive *Listeria* colonies are selected under 45° transillumination. *Listeria* colonies are dense white to iridescent white appearing as crushed glass.

Small colonies tend to be blue, while non-*Listeria* show yellowish orange colonies that are further purified.

## Methodology

Suspend 46.00 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool below 50°C. Before gelling, aseptically add sterile rehydrated contents of 1 vial of McBride Listeria Supplement (MS2070). Shake well before pour into sterile Petri plates.

**Warning:** Lithium chloride is harmful. Avoid bodily contact and inhalation of vapours. On contact with skin, wash with plenty of water immediately.

## Quality Control

### Appearance

Cream to yellow homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel

### Colour and Clarity

Light amber clear to slightly opalescent gel forms in Petri plates

### Reaction

Reaction of 4.6% w/v aqueous solution at 25°C. pH : 7.3±0.2

### pH Range

7.10-7.50

### Cultural Response

DM1891: Cultural characteristics observed under anaerobic condition with added McBride Listeria Supplement (MS2070) after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Enterococcus faecalis</i> ATCC 29212	50-100	none – poor	0-10%
<i>Escherichia coli</i> ATCC 25922	50-100	good-luxuriant	>=50%
<i>Listeria monocytogenes</i> ATCC 19112	50-100	none – poor	0-10%
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	none – poor	0-10%

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Hyslop N., St. G. and Osborne A. D., 1959, Vet. Rec. 71 : 1082
2. Gray M. L., and McWade D. H., 1964, J. Bacteriol., 68:634
3. Murray E. G. D., Webb R. A. and Swann M. B. R., 1926, J. Pathol. Bacteriol., 29:407.
4. Conner D. E., Brackett R. E., and Beuchat L. R., 1986, Appl. Environ. Microbiol., 52: 59
5. Doyle M. P., Meske L. M. and Marth E. H., 1985, J. Food Prot., 48:740
6. Dowell, Hill and Altemeier, 1964, J. Bacteriol., 88:1811.
7. FDA Bacteriological Analytical Manual, 2005, 18th Ed., AOAC, Washington, DC.
8. McClain D. and Lec W., 1989, Lab. Comm No. 57, Revised, May 24, 1989, U.S. Dept. of Agric., FSIS, Microbiol. Div., Beltsville, MD.



Dehydrated Culture Media  
Bases / Media Supplements

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performances parameters.

