

## Technical Information

### Kundrat Agar

#### Product Code: DM 2360

**Application:** - Kundrat Agar is recommended for the routine qualitative detection of residues from antibiotics and other chemotherapeutical agents in animal-derived food

#### Composition\*\*

| Ingredients         | Gms / Litre |
|---------------------|-------------|
| Meat peptone        | 7.800       |
| Casein peptone      | 7.800       |
| Yeast extract       | 2.800       |
| Sodium chloride     | 3.000       |
| Dextrose            | 1.000       |
| Starch              | 4.000       |
| Gelatin             | 4.000       |
| Bromocresol purple  | 0.016       |
| Agar                | 10.000      |
| Final pH ( at 25°C) | 6.8±0.2     |

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Kundrat Agar is used for detection of antimicrobial residues in animal feed preparations. The test is carried out using a Ampoule of *Bacillus stearothermophilus* (LA840) as test microorganisms. It is also used for detection of antimicrobial residues in meat and organ samples; used together with spore suspensions of *Bacillus subtilis* (BGA) as test organism. Presence of chemotherapeutic agents is indicated by the formation of inhibition halos or zones around the disc with the sample(1).

The test is performed in the form of an agar diffusion test. Any inhibitors present produce inhibition zones devoid of bacterial growth surrounding the applied samples. With further incubation, the test organism ferments glucose present in the medium to form acid, that causes bromocresol purple to change its colour to yellow. Only the inhibition zone still retains the original violet colour of the indicators.

**Test Procedure:** After autoclaving the medium, cool to 50-60°C. To each 200 ml of the medium add the contents of 1 ampoule of *B.stearothermophilus* Ampoule (LA840), mix, pour in plates.

Filter paper discs with a diameter of 6 mm are soaked with the liquid specimen or placed on organ (kidney, liver) or muscle sections. The discs are then slightly pressed onto the surface of the culture medium (up to 6 discs per plate) (2).

Two methods are recommended for performing the test:

1)45 minutes incubation, rapid test:

After placing the discs on the preincubated plates, incubate them for further 45 minutes at 65°C without prediffusion.

2)3 hour incubation:

The plates are not preincubated. After the filter paper discs have been applied to the plates, they should be incubated for 3 hours at 65°C without pre-diffusion.

In the case of rapid test, formation of inhibition zones can be seen after 15-25 minutes incubation in the medium, which is otherwise turbid as a result of spore growth. After the 45 minutes incubation, the inhibition zones become even more distinct due to the fact that the culture medium changes colour. Formation of inhibition zones is to be regarded as a positive result.



Dehydrated Culture Media  
Bases / Media Supplements

In the case of the 3 hours incubation, only those inhibition zones with a diameter of more than 10 mm can be considered positive. If a distinct colour change has not occurred after 45 minutes or 3 hours, incubation can be prolonged.

Cleaning agents, disinfectants and preservatives are not covered by this test. When performing the rapid test, pre-incubating the inoculated plates enhances growth of the test organism; the inhibition zones then appear more rapidly after application of the samples.

## Methodology

Suspend 40.41 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat to boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well before pouring into sterile Petri plates.

## Quality Control

### Appearance

Cream to yellow coloured with green tinge, homogeneous free flowing powder

### Gelling

Firm, comparable with 1.0% Agar gel.

### Colour and Clarity

Light purple coloured, clear to slightly opalescent gel forms in Petri plates

### Reaction

Reaction of 4.04% w/v aqueous solution at 25°C. pH : 6.8±0.2

### pH Range

6.60-7.00

### Cultural Response

DM2360: Cultural characteristics observed after an incubation at 65°C for 18-24 hours.

### Organism

### Growth

### Cultural Response

*Bacillus stearothermophilus* ATCC 7953

good-luxuriant

## Storage and Shelf Life

**Dried Media:** Store below 30°C and prepared medium at 2-8°C. Use before expiry period on the label

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Kundrat W., 1968, Methoden zur Bestimmung von Antibiotika-Rückständen in tierischen Produkten. - Z. Anal. Chem.; 624-630.
2. Kundrat W., 1972, 45- Minuten - Schnellmethode zum mikrobiologischen Nachweis von Hemmstoffen in tierischen Produkten. - Fleischwirtsch., 52; 485-487.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performances parameters.

