

Technical Information

Clostridium Broth Base

Product Code: DM 2315

Application: - Clostridium Broth Base is used for identification of spores of *Clostridium tyrobutyricum*, which is usually responsible for late blowing in cheese.

Composition**

Ingredients	Gms / Litre
Casein peptone	15.000
Meat extract	10.000
Yeast extract	5.000
Sodium acetate	5.000
L-Cysteine	0.500
Final pH (at 25°C)	6.0±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Ripening of cheese under controlled conditions of temperature and humidity determine the final flavour and body characteristics of the product. Microbial spoilage in cheese is generally limited because of the combined effect of acid and salt and is less likely in the lower moisture cheese. Spores of *Clostridium tyrobutyricum* in the milk used for the manufacture of Emmentaler, Edam and Gouda can survive the heat treatment used for cheese milk and cause late gas formation (blowing defect) and related off-flavours during ripening (1). Even low spore densities of this anaerobe in milk used for cheese production can bring about this phenomenon, if the growth conditions are suitable. Clostridium Broth Base is used for the identification of spores of *C. tyrobutyricum*.

C. tyrobutyricum ferments lactate and acetate to butyrate, CO₂ and H₂ (2, 3). As against *C. butyricum*, *C. tyrobutyricum* grows in media with lower acidic pH (4, 5). Low pH value of the media helps to inhibit the other microbial flora thereby favoring growth of *C. tyrobutyricum*.

Casein peptone, meat extract and yeast extract supply the essential nutrients mainly the nitrogen compounds whereas L-cysteine promotes the growth of Clostridial species by maintaining a low oxygen tension in the medium.

Methodology

Suspend 35.5 grams of dehydrated media in 1000 ml distilled water. Mix thoroughly & heat if necessary to dissolve the medium completely. Add 10 ml of 50% sodium lactate. Mix well and dispense into tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity

Yellowish brown coloured, clear solution without any precipitate

Reaction

Reaction of 3.55% w/v aqueous solution at 25°C. pH : 6.0±0.2

Ph Range

5.80-6.20

Cultural Response

DM 2315: Cultural characteristics observed under anaerobic condition, after an incubation at 35-37°C for up to 7 days.



Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth	Gas
Cultural Response			
<i>Clostridium perfringens</i> ATCC 10543	50-100	luxuriant	positive reaction
<i>Clostridium tyrobutyricum</i> ATCC 25755	50-100	luxuriant	positive reaction
<i>Escherichia coli</i> ATCC 25922	50-100	good	variable reaction
<i>Pseudomonas aeruginosa</i> ATCC 27853	$\geq 10^3$	inhibited	-
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good	negative reaction

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Vanderzant C. and Splittstoesser D. F., (Eds.), 1992, Compendium of Methods of for the Microbiological Examination of Foods, 3rd Ed., APHA, Washington, D.C.
2. Goudkov A. V. and Sharp, M. E., 1966, J. Dairy Res., 33: 139-149.
3. Bryant M. P. and Burkey L. A., 1956, J. Bacteriol., 71: 43-46.
4. Kutzner H. J., 1963, Infektionskr. Hyg. Abt. 1, 191: 441-450.
5. Balows A. et al, 1992, The Prokaryotes, A Handbook of the Biology of Bacteria : Ecophysiology, Isolation , Identification Applications, 2nd Ed. Vol. II. p. 1800-1866.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performances parameters.

