

Technical Information

Letheen Broth w/Triton X-100

Product Code: DM 2459

Application: - Letheen Broth w/Triton X - 100 is used for screening cosmetic products for microbial contamination.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Beef extract	5.000
Sodium chloride	5.000
Lecithin	0.700
Polysorbate 80	5.000
Triton X-100	1.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Letheen Broth with Triton X-100 is recommended for luxuriant growth of most organisms for detection of yeast and moulds. Triton X-100 is non-ionic and disperses microorganisms making counting easier.

In the early 40s, Weber and Black recommended the use of lecithin and polysorbates to neutralize the antimicrobial action of the quaternary ammonium compounds (4). In 1965, the methodology was accepted by AOAC for the antimicrobial assays and extended their use to all the cationic detergents. In 1978, the FDA incorporated it as pre-enrichment medium for every microbial examination of cosmetics.

There are great chances of altering the chemical composition of cosmetics by the metabolism of organisms thereby spoiling and causing harm to the users (1, 2, 3). Direct colony counts and enrichment culturing are the methods of choice for isolating microorganisms from cosmetic products. The word Letheen represents a combination of lecithin and polysorbate (tween) 80.

Peptic digest of animal tissue, beef extract supply nitrogenous nutrients, carbon compounds and trace elements for the growth of microorganisms. Incorporation of lecithin and polysorbate 80 to the medium enables the recovery of bacteria from materials containing residues of disinfectant compounds or preservatives used in cosmetics. Polysorbate 80 is added to nullify phenolic compounds, hexachlorophene, formalin and along with lecithin neutralizes ethyl alcohol (5). Lecithin also neutralizes quaternary ammonium compounds present in the cosmetics. Sodium chloride maintains the osmotic balance of the medium. Triton X-100 acts as a surfactant. Cosmetics contain preservatives and they should be at least partially inactivated during the plating and this medium helps in dilution as well as neutralizing.

Methodology

Suspend 26.7 grams of dehydrated powder media in 1000 ml distilled water. Mix thoroughly & heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Shake well before dispense as desired.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder



Dehydrated Culture Media
Bases / Media Supplements

Colour and Clarity

Yellow coloured, clear solution in tubes

Reaction

Reaction of 2.67% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH Range

6.80-7.20

Cultural Response

DM 2459: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant
<i>Staphylococcus aureus</i> ATCC 6538	50-100	luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Dunningan A. P., 1968, Drug Cosmet. Ind., 102:43.
2. Smart R. and Spooner D. F., 1972, J. Soc. Cosmet. Chem., 23:721.
3. Wilson L. A. and Ahearn D. G., 1977, Am. J. Ophthalmol., 84:112.
4. Weber and Black, 1948, Soap Sanitary Chem., 24:134-139
5. Favero (Chm.), 1967, A State of the Art Report, Biological Contamination Control Committee, American Association for Control.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

