



ANTIMICROBIAL SUSCEPTIBILITY
SYSTEMS

Technical Information

Sterile Discs*

Product Code: SD1067

Application: Sterile discs can be used for testing antibiotic susceptibility testing by loading antibiotic of choice.

Appearance: Sterile filter paper discs of 6mm diameter.

Cultural response: Sterility of the discs checked for 14 days in Alternative Thioglycollate Medium (DM1010) at 30-35°C and Soyabean Casein Digest Medium (DM1011) at 20-25°C.

Sample	Results in Alternative Thioglycollate Medium (DM010)	Results in Soyabean Casein Digest Medium (DM1011)
Test	Clear, No growth observed	Clear, No growth observed
Negative Control	Clear, No growth observed	Clear, No growth observed
Positive control	Luxuriant growth observed	Luxuriant growth observed

Storage and Shelf Life

Store at 2°C to 8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.