



# Product Specification

cdhfinechemical.com

## CEDEPOL™ ALSO KNOWN AS TEEPOL

<b>PRODUCT CODE</b>	<b>056011</b>
<b>SYNONYMS</b>	--
<b>C.I. NO.</b>	--
<b>CASR NO.</b>	--
<b>ATOMIC OR MOLECULAR FORMULA</b>	--
<b>ATOMIC OR MOLECULAR WEIGHT</b>	--
<b>PROPERTIES</b>	--

PARAMETER	LIMIT
Description	Pale yellow, almost clear liquid.
Solubility	5% Solution in water is almost clear and colourless.
pH of (10% solution in water)	6.0-7.5
wt. per ml at 20 <sup>o</sup> C	About 1.06 g
Suitability test	Passes test

(Use as cleansing and degreasing agent in Laboratories, Hospitals, Pathological diagnostic centers etc. For cleaning Glass and porcelain apparatus, plastic and metals containers, slides lenes Tiles, Working benches, Floor and walls )Use by mixing 5-10 mk of cedepol with 1 litre of water. For Lenses. Wood work and other delicate work ise 5-10 ml cedepol with 1 litre of cold water.

**Note(s) : Assay (if applicable) method mentioned.**

<p><b>WARNING</b></p> <p><b>Hazard statements:-</b> Not hazardous. No hazards.</p> <p><b>Precautionary statements</b></p> <p><b>Prevention:</b> -----</p> <p><b>Response:</b> -----</p> <p><b>Disposal:</b>-----</p> <p><b>Hazard Pictogram(s) :</b>-----</p>	<p>IMDG Code : --</p> <p>UN No. : --</p> <p>IATA : --</p>
---	---