



## **Product Specification**

cdhfinechemical.com

	CAPSULE STAIN KIT
PRODUCT CODE	K010005
SYNONYMS	
C.I. NO.	
CASR NO.	<del></del>
ATOMIC OR MOLECULAR FORMULA	
ATOMIC OR MOLECULAR WEIGHT	
PROPERTIES	<del></del>
COMPOSITION	
Methylene Blue (aqueous) (865650)	
Ingredients	
Methylene blue	0.3 gm
Distilled water	100.0ml
Nigrosin Stain 10%w/v (869200)	
Ingredients	
Nigrosin	10.0 gm
Distilled water	100.0 ml
Formalin	0.5 ml

WARNING

IMDG Code : 9/II UN No. HAZARD STATEMENTS: May cause respiratory irritation. Flammable liquid and vapour. Harmful if swallowed. Causes mild : 3316 skin irritation. Causes eye irritation. IATA : 9

## PRECAUTIONARY STATEMENTS:

Prevention: Ground/bond container and receiving equipment. Use only nonsparking tools. Take precautionary Measures against static discharge. Wash hands thoroughly after handling. Do not eat, drink or smoke when using this product. Keep away from heat/sparks/open flame - No smoking. Keep container tightly closed. Use explosion-proof electrical/ventilating/lighting/equipment. Wear protective gloves and eye/face protection.

Response: If eye irritation persists, get medical advice/attention. If skin irritation occurs, seek medical advice/attention. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment: refer to Label or MSDS. If on skin or hair: emove/take off immediately all contaminated clothing. Rinse with water/shower. In case of fire, use foam for extinction.

Disposal: Allow to evaporate in safe open area (applicable for small quantities). The chemical should be burnt in a chemical ncinerator equipped with burner and scrubber. If the volumes are larger than 500ml and no incinerator and site of evaporation is not available then transfer the waste chemical to waste solvent drum for collection by specialist disposal company.

Hazard Pictogram(s):





Replace Date 20-June-2023